

Complicated Cataracts, Simplified



READY-LOADED IRIS SPECULUM



Xpand^{NT} Insertion


Click play above to watch the X1 Ready-Loaded Usage Guide



DIAMATRIX[®]

What is that awesome device?

The X1 system is unlike any iris expansion product on the market. Ready-Loaded with the XpandNT Iris Speculum, the X1 was developed expert in conjunction with industry leaders including surgeons and administrators. Among the many realized benefits of this device is that the device is atraumatic to ocular structures and it is easy to employ. Check out the other key benefits.



Key Benefits

1. Nitinol Ring

6.7mm aperture created from a laser-welded wire. Smooth and Flexible; Biocompatible; .003" diameter wire.

2. Low Profile

.8mm pockets gently cradle the iris rim.

3. 8 points of contact


More contact provides less iris trauma and a more circular pupil.

4. Speculum design

Creates stability and enables easy rotation.

How do I get that awesome device in the eye?

With our innovative injector! After years of customer engagement, we developed an injector that works so well with the XpandNT, you might think about keeping it. Alas, the X1 is a single-use



Still Have Questions?

We know people. If you have questions, please call us at 800.867.8081 and speak with our customer service agent. Don't have time to talk, email us cs@diamatrix.com