



DIAMATRIX

Ophthalmic Instruments

A LEGACY OF
Quality, Value, and Innovation



WARRANTY TERM

The legacy of producing innovative quality instruments continues today at Diamatrix Ltd. We unconditionally guarantee our ophthalmic instruments to be of the highest quality in craftsmanship and free of material and manufacturing defects for the lifetime of the instrument.

WARRANTY COVERAGE

Diamatrix Ltd will provide lifetime replacement or repair of our instruments, provided that the instrument is used for its intended surgical purpose. The expected life of each instrument varies by type. We will either repair or replace defective instruments without charge at our discretion.

Any intentional or unintentional changes, alterations, modifications, or non-Diamatrix repairs will release us from this warranty. Diamatrix Ltd instrument warranty is non-transferable.

The Diamatrix Instrument Lifetime Warranty only applies to Diamatrix Instruments and does not extend to our knives, suture, or other products. A partial list of included instruments is included below:

- Forceps
- Scissors
- Handheld Instruments (item number beginning with "H")
- I/A and Bimanual Handpieces
- Markers
- Needle Holders

LIMITATIONS AND RESTRICTIONS

The following limitations and restrictions apply to this warranty:

- 1.** The Warranty includes only those products sold as new from Diamatrix Ltd.
- 2.** Indications that the instrument was not used for its intended surgical purpose or was not cleaned in a manner consistent with our Directions of Use shall void the warranty.

3. Repairs, modifications, or alternations to the instrument made by any other entity, will breach this lifetime warranty.

4. This guarantee does not apply to instruments that are mishandled, misused, improperly maintained, improperly sharpened, or improperly sterilized.

Please refer to our product info tab on our website for complete details of proper handling and care of Diamatrix Ophthalmic Instruments.

REPAIRS OF NON-DIAMATRIX INSTRUMENTS

Diamatrix does offer a full repair service for non-Diamatrix branded instruments, or for those instruments that are no longer within their "useful life".

We do offer a 1-year warranty on all repair work and ensure that your repaired instrument is in a like new condition.

Upon inspection our highly skilled technicians will determine if your instrument can be repaired to function comparable to a new instrument or if the instrument should be replaced. Instruments that have been repaired or sharpened by a source other than Diamatrix Ltd are void of this warranty.

RETURN POLICY

If you are not not 100% satisfied with the instrument, contact Diamatrix within 15 calendar days of physical delivery date and Diamatrix will issue you a full refund. Product must be returned without damage due to mishandling or improper use. A credit will be issued to the customer account only after the item to be returned has been received. If items are not returned to Diamatrix, payment will still be due.

To initiate a return, please contact Diamatrix Customer Service at 1-800-867-8081. A Return Goods Authorization (RGA) number will be issued. Please reference your RGA number on the outside of your shipping container and send to Diamatrix at:

Diamatrix
210 Nursery Road
The Woodlands, TX 77380

CARE, CLEANING, AND STERILIZATION INSTRUCTIONS:

CAUTION: Incorrect care, cleaning, and misuse can damage delicate surgical instruments. Damage to the instrument may lead to patient harm. Please read these instructions prior to use. Users should be trained in the use and handling of surgical instruments.

WARNING: EVERY INSTRUMENT MUST BE CLEANED AND STERILIZED PRIOR TO FIRST USE.

Instructions prior to use:

Users must inspect the instrument prior to use or reprocessing for breaks, cracks, malformations of material, or malfunctions. Do not use damaged instruments. Do not attempt to repair any damaged instruments. Repairs should be made by Diamatrix and its subcontractors only. Questions about repairs of Diamatrix instruments should be forwarded to Diamatrix Customer Service.

Cleaning:

Immediately after surgery, rinse the instrument with distilled water. Never let foreign materials, such as BSS or proteins, dry on the instrument. Note: Based on the documented risk of TASS associated with enzyme detergent use, without proven benefit for endophthalmitis prevention, enzymatic detergent should not be required for routine decontamination of ophthalmic intraocular instruments.

Clean instruments with mild soap and distilled water. Cleaning solutions and rinses at or near a neutral pH (7.0) are best. Do not use abrasive cleaners.

Flush thoroughly with distilled/deionized water.

Dry instruments completely and carefully with compressed air (including inside lumen and other areas that may hold moisture).

Ultrasonic Cleaning:

The best method for cleaning the instrument is to use an ultrasonic cleaner. Ultrasonic cleaning is recommended for all instruments at the end of each surgical day or after a number of surgical procedures as indicated by the appearance of the instrument.

Follow the instructions as supplied by the manufacturer of your ultrasonic cleaner. Always use distilled water for ultrasonic cleaning or a cleaning agent specifically recommended for surgical instruments.

A water temperature of approximately 150°F (65°C) is ideal for this purpose and will produce the best cleaning results. Many ultrasonic cleaners come with built-in heaters, but if your unit does not have one, you may use an immersion heater or heat the water in a separate stainless-steel container.

The instrument(s) should be placed in a wire or perforated plastic basket and suspended in the cleaning solution to avoid damaging delicate tips. Instruments must be totally submerged in the solution but should not be in contact with each other. If you use ultrasonic cleaning regularly, a cleaning cycle of approximately 5 minutes should be adequate, but some debris may require additional time.

After the ultrasonic cycle, thoroughly rinse the instrument under running water and then follow with a final rinse in a clean bath of distilled water. Cleaning and rinsing solutions should be replaced frequently.

Air-dry the instrument or dry with a hot-air blower or a lint-free cloth before returning it to storage.

Most instruments should not require lubricating baths, but it is recommended that an instrument with joints and moving parts be lubricated after cleaning. Treat such parts with a steam-permeable, water-soluble, instrument lubricant following each cleaning and prior to sterilization. Sterilization **CAUTION:** Autoclave temperatures should not exceed 280°F (137°C) to avoid damage to any nonmetallic parts.

Note: If applicable, instruments should be sterilized in the open or unlocked position.

Note: Instruments made of different alloys should be cleaned and sterilized separately.

Steam autoclave sterilization is recommended.

Standard autoclave cycle:

Steam sterilize at 270°F (132°C) for 5 minutes.

Note: Other time and temperature cycles may be used but the user must validate any deviation from the recommended time and temperature. 20 minutes of dry time is determined to be sufficient to completely dry the instruments. Refer to the manufacturer of your autoclave to confirm appropriate cycle settings.

Inspection and Storage:

All microsurgical instruments should be inspected under magnification at the end of each surgical day by qualified personnel to ensure they're in proper working condition. Damaged instruments should be removed from circulation and immediately repaired or replaced. Instruments with delicate tips should be stored using a protective cover.

The instrument should be stored in the same container in which it will be sterilized. If multiple instruments are in the same tray, they should not be touching.

Sales Support & Representation:

- Diamatrix will cross-reference to any competitor's part number and provide quotations.
- For a Sales Representative or Diamatrix Distributor in your area, please call 800-867-8081

Repairs:

- Diamatrix provides a full instrument repair service. To find out more information, please call: 800-867-8081
- Please call for a Returned Good Authorization Number.
- Instruments and diamond knives have a one (1) year warranty against defects in material and workmanship, which applies from the date received.
- Damage resulting from misuse or mishandling, modifications made by another party void Diamatrix warranty.








Bimanual Hand Pieces and Tips.....	1
I&A Tips Hand Pieces and Tips.....	12
Cannulas.....	20
Fixation Rings.....	70
Forceps.....	74
Handheld Instruments.....	109
Markers.....	185
Needle Holders.....	212
Scissors.....	218
Wire Speculum.....	234



Bimanual Hand Pieces and Tips





Bimanual Hand Pieces and Tips

B-00

O-Rings

Bimanual Replacement O-Rings (5 Pieces)



B1-00

Diamatrix

Bimanual Replacement Blue Irrigating Handle Titanium



B1-21-02

Diamatrix


Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Curved shaft, .5mm dual ports, Round Irrigating blue handle, 111mm total length



B1-21-02T

Diamatrix


Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Curved shaft, .5mm dual ports, 21mm tip length



B1-21-03

Diamatrix

Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Curved shaft, .4mm dual ports, Round Irrigating blue handle, 111mm total length



Bimanual Hand Pieces and Tips

B1-21-03T

Diamatrix


Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Curved shaft, .4mm dual ports, 21mm tip length



B1-21-04

Diamatrix


Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Curved shaft, .3mm dual ports, Round Irrigating blue handle, 111mm total length



B1-21-04T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Curved shaft, .3mm dual ports, 21mm tip length



B1-21-05

Diamatrix

Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Straight shaft, Open front, Round Irrigating blue handle, 111mm total length



B1-21-05T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Straight shaft, Open front, 21mm tip length



Bimanual Hand Pieces and Tips

B1-21-06

Diamatrix


Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Curved shaft, Open front, Round Irrigating blue handle, 111mm total length



B1-21-06T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Curved shaft, Open front, 21mm tip length



B1-21-07

Diamatrix

Bimanual Irrigating Blue Handpiece, 21 Gauge Cannula tip, Straight shaft, .5mm dual ports, Round Irrigating blue handle, 111mm total length



B1-21-07T

Diamatrix


Bimanual Irrigating Cannula Replacement tip, 21 Gauge, Straight shaft, .5mm dual ports, 21mm tip length



B1-22-15

Diamatrix

Bimanual Irrigating Blue Handpiece, 22 Gauge Cannula tip, Curved shaft, Open front, Round Irrigating blue handle, 111mm total length



Bimanual Hand Pieces and Tips

B1-22-15T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 22 Gauge, Curved shaft, Open front, 21mm tip length



B1-23-01

Diamatrix

Bimanual Irrigating Blue Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm dual ports, Round Irrigating blue handle, 111mm total length



B1-23-01T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm dual ports, 21mm tip length



B1-23-03

Diamatrix


Bimanual Irrigating Blue Handpiece, 23 Gauge Thin Wall Cannula tip, Straight shaft, .4mm dual ports, Round Irrigating blue handle, 111mm total length



B1-23-03T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 23 Gauge Thin Wall, Straight shaft, .4mm dual ports, 21mm tip length



Bimanual Hand Pieces and Tips

B1-23-04

Diamatrix

Bimanual Irrigating Blue Handpiece, 23 Gauge Thin Wall Cannula tip, Curved shaft, .4mm dual ports, Round Irrigating blue handle, 111mm total length



B1-23-04T

Diamatrix


Bimanual Irrigating Cannula Replacement tip, 23 Gauge Thin Wall, Curved shaft, .4mm dual ports, 21mm tip length



B1-23-08

Diamatrix


Bimanual Irrigating Blue Handpiece, 23 Gauge Thin Wall Cannula tip, Curved shaft, .4mm dual oval ports, Round Irrigating blue handle, 111mm total length



B1-23-08T

Diamatrix


Bimanual Irrigating Cannula Replacement tip, 23 Gauge Thin Wall, Curved shaft, .4mm dual oval ports, 21mm tip length



B1-23-12

Diamatrix

Bimanual Irrigating Blue Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm single port, Round Irrigating blue handle, 111mm total length



Bimanual Hand Pieces and Tips

B1-23-12T

Diamatrix

Bimanual Irrigating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm single port, 21mm tip length



B2-00

Diamatrix


Bimanual Replacement Gold Aspirating Handle Titanium



B2-21-03

Diamatrix

Bimanual Aspirating Gold Handpiece, 21 Gauge Cannula tip, Curved shaft, .3mm single port, Round Aspirating gold handle, 111mm total length



B2-21-03T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 21 Gauge, Curved shaft, .3mm single port, 21mm tip length



B2-21-04

Diamatrix

Bimanual Aspirating Gold Handpiece, 21 Gauge Cannula tip, Curved shaft, .3mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length




Bimanual Hand Pieces and Tips

B2-21-04T

Diamatrix

Bimanual Aspirating Cannula Replacement tip, 21 Gauge, Curved shaft, .3mm single port with Capsule Polisher tip, 21mm tip length



B2-21-05

Diamatrix


Bimanual Aspirating Gold Handpiece, 21 Gauge Cannula tip, Curved shaft, .4mm single port, Round Aspirating gold handle, 111mm total length



B2-21-05T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 21 Gauge, Curved shaft, .4mm single port, 21mm tip length



B2-21-06

Diamatrix


Bimanual Aspirating Gold Handpiece, 21 Gauge Cannula tip, Curved shaft, .4mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



B2-21-06T

Diamatrix

Bimanual Aspirating Cannula Replacement tip, 21 Gauge, Curved shaft, .4mm single port with Capsule Polisher tip, 21mm tip length




Bimanual Hand Pieces and Tips

B2-21-08

Diamatrix


Bimanual Aspirating Gold Handpiece, 21 Gauge Cannula tip, Curved shaft, .35mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



B2-21-08T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 21 Gauge, Curved shaft, .35mm single port with Capsule Polisher tip, 21mm tip length



B2-22-10

Diamatrix


Bimanual Aspirating Gold Handpiece, 22 Gauge Cannula tip, Curved shaft, .35mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



B2-22-10T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 22 Gauge, Curved shaft, .35mm single port with Capsule Polisher tip, 21mm tip length



B2-23-01

Diamatrix

Bimanual Aspirating Gold Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm single port, Round Aspirating gold handle, 111mm total length




Bimanual Hand Pieces and Tips

B2-23-01T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm single port, 21mm tip length



B2-23-02

Diamatrix


Bimanual Aspirating Gold Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



B2-23-02T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm single port with Capsule Polisher tip, 21mm tip length



B2-23-03

Diamatrix


Bimanual Aspirating Gold Handpiece, 23 Gauge Thin Wall Cannula tip, Curved shaft, .4mm single port with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



B2-23-03T

Diamatrix

Bimanual Aspirating Cannula Replacement tip, 23 Gauge, TW, Curved shaft, .4mm single port with Capsule Polisher tip, 21mm tip length




Bimanual Hand Pieces and Tips

B2-23-04

Diamatrix

Bimanual Aspirating Gold Handpiece, 23 Gauge Thin Wall Cannula tip, Curved shaft, .4mm single port, Round Aspirating gold handle, 111mm total length



B2-23-04T

Diamatrix


Bimanual Aspirating Cannula Replacement tip, 23 Gauge, TW, Curved shaft, .4mm single port, 21mm tip length



B2-23-07

Diamatrix

Bimanual Aspirating Gold Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm dual oval ports, Round Aspirating gold handle, 111mm total length



B2-23-07T

Diamatrix

Bimanual Aspirating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm dual oval ports, 21mm tip length



B2-23-13

Diamatrix

Bimanual Aspirating Gold Handpiece, 23 Gauge Cannula tip, Curved shaft, .3mm dual oval ports with Capsule Polisher tip, Round Aspirating gold handle, 111mm total length



Bimanual Hand Pieces and Tips

B2-23-13T

Diamatrix

Bimanual Aspirating Cannula Replacement tip, 23 Gauge, Curved shaft, .3mm dual oval ports with Capsule Polisher tip, 21mm tip length



B2-23-14

Randolph

Bimanual Aspirating Gold Handpiece, 23 Gauge Randolph Cannula tip, Spatulated shaft, Angled 12mm from bend to tip, Round Aspirating gold handle, 111mm total length



B2-23-14T

Randolph

Bimanual Aspirating Cannula Replacement tip, 23 Gauge, Randolph Cannula tip, Spatulated shaft, Angled 12mm from bend to tip, 25mm tip length





I&A Hand Pieces and Tips



I&A Hand Pieces and Tips

IS-ASP

Diamatrix


Aspirating Luer Adaptor for easy use with any Standard Aspiration Cannula, Used on the Irrigating/Aspirating Handpieces, Stainless



IS-H

Diamatrix

Diamatrix Irrigating and Aspirating Handpiece Handle, 105mm total length, Diamatrix Signature Feather-Series round handle, Titanium for Stainless style tips



I-SIL

Diamatrix


Replacement Silicone Sleeves for Irrigating/Aspirating Handpiece tips Pack of 10 pieces



IS-IRR

Diamatrix

Irrigating Luer Adaptor for easy use with any Standard Irrigation Cannula, Used on the Irrigating/Aspirating Handpieces, Stainless



IS-O

Diamatrix

Replacement O-ring sets for Stainless Handpieces
Pack contains 5 of each



5 pieces of each O'ring

IS-3-C

Diamatrix

Irrigating/Aspirating Curved tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 16mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3-S

Diamatrix

Irrigating/Aspirating Straight tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3S-C

Diamatrix

Irrigating/Aspirating Silicone Curved tip only, Silicone Outer Sleeve, 19 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 16mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3S-S

Diamatrix

Irrigating/Aspirating Silicone Straight tip only, Silicone Outer Sleeve, 19 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainl



IS-3S-45

Diamatrix

Irrigating/Aspirating Silicone 45° tip only, Silicone Outer Sleeve, 19 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3-20

Diamatrix

Irrigating/Aspirating 20° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3-45

Diamatrix

Irrigating/Aspirating 45° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3-90

Diamatrix

Irrigating/Aspirating 90° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-3-150

Diamatrix

Irrigating/Aspirating 150° tip only, 16 Gauge Outer Sleeve, 23 Gauge Inner Cannula tip with .3mm front port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-5-C

Diamatrix

Irrigating/Aspirating Curved tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 16mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-5-S

Diamatrix

Irrigating/Aspirating Straight tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-5S-45

Diamatrix

Irrigating/Aspirating Silicone 45° tip only, Silicone Outer Sleeve, 19 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-5-20

Diamatrix

Irrigating/Aspirating 20° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless

No Image Available

IS-5-45

Diamatrix

Irrigating/Aspirating 45° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IS-5-90

Diamatrix

Irrigating/Aspirating 90° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



I&A Hand Pieces and Tips

IS-5-150

Diamatrix


Irrigating/Aspirating 150° tip only, 16 Gauge Outer Sleeve, 21 Gauge Inner Cannula tip with .5mm front port, 17mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Stainless



IT-ASP

Diamatrix

Aspirating GOLD Luer Adaptor for easy use with any Standard Aspiration Cannula, Used on the Irrigating/Aspirating Handpieces, Titanium



IT-H

Diamatrix


Diamatrix Irrigating and Aspirating Handpiece Handle, 105mm total length, Diamatrix Signature Feather-Series round handle, Titanium



IT-IRR

Diamatrix

Irrigating BLUE Luer Adaptor for easy use with any Standard Irrigation Cannula, Used on the Irrigating/Aspirating Handpieces, Titanium



IT-O

Diamatrix

Replacement O-ring sets for Titanium Handpieces Pack contains 5 of each



5 pieces of each O'ring

I&A Hand Pieces and Tips

I-TEST

Diamatrix

Replacement Test Chambers for Irrigating/Aspirating Handpieces Pack of 5 pieces



IT-3-C

Diamatrix

Irrigating/Aspirating Curved tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 18mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3-S

Diamatrix

Irrigating/Aspirating Straight tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3S-C

Diamatrix


Irrigating/Aspirating Silicone Curved tip only, Silicone Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 18mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3S-S

Diamatrix

Irrigating/Aspirating Silicone Straight tip only, Silicone Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3S-45

Diamatrix

Irrigating/Aspirating Silicone 45° tip only, Silicone Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3-45

Diamatrix

Irrigating/Aspirating 45° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3-90

Diamatrix

Irrigating/Aspirating 90° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .3mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-3-150

Diamatrix

Irrigating/Aspirating 150° tip only, 16 Gauge Outer Sleeve, 23 Gauge Inner Cannula tip with .3mm front port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-4-S

Diamatrix

Irrigating/Aspirating Straight tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .4mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-5-C

Diamatrix

Irrigating/Aspirating Curved tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 18mm Curved Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-5-S

Diamatrix

Irrigating/Aspirating Straight tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-5-45

Diamatrix

Irrigating/Aspirating 45° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-5-90

Diamatrix

Irrigating/Aspirating 90° tip only, 16 Gauge Outer Sleeve, 21 Gauge Closed Inner Cannula tip with .5mm Aspirating Port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing



IT-5-150

Diamatrix

Irrigating/Aspirating 150° tip only, 16 Gauge Outer Sleeve, 21 Gauge Inner Cannula tip with .5mm front port, 19mm Straight Shaft, Female Connector on Irrigation Line and Male Connector on Aspirating Line, 37mm total length, Titanium with Stainless Inner Tubing





Cannulas




Cannulas

C027-8

Diamatrix


Air Injection 27 Gauge, Angled 8mm from bend to tip,
19mm tip length



C030-5

Diamatrix


Air Injection 30 Gauge, Angled 5mm from bend to tip,
19mm tip length



C030-6

Diamatrix


Air Injection 30 Gauge, Angled 6mm from bend to tip,
19mm tip length



C030-7

Diamatrix


Air Injection 30 Gauge, Angled 7mm from bend to tip,
19mm tip length



C030-8

Diamatrix

Air Injection 30 Gauge, Angled 8mm from bend to tip,
19mm tip length




Cannulas

C119-19

Viscoelastic


Viscoelastic Aspirating Cannula, 19 Gauge, Angled
10mm from bend to tip, 25mm tip length



C119-21

Viscoelastic


Viscoelastic Aspirating Cannula, 21 Gauge, Angled
7mm from bend to tip, 25mm tip length



C119-22

Viscoelastic


Viscoelastic Aspirating Cannula, 22 Gauge, Angled
10mm from bend to tip, 25mm tip length



C119-26

Viscoelastic


Viscoelastic Aspirating Cannula, 26 Gauge, Angled
10mm from bend to tip, 25mm tip length



C125-23

Rycroft

Rycroft Anterior Chamber Cannula, 23 Gauge, Angled
45° 5mm from bend to tip, 22mm tip length




Cannulas

C125-25

Rycroft


Rycroft Anterior Chamber Cannula, 25 Gauge, Angled 45° 5mm from bend to tip, 22mm tip length



C125-27

Rycroft


Rycroft Anterior Chamber Cannula, 27 Gauge, Angled 45° 5mm from bend to tip, 22mm tip length



C125-30

Rycroft


Rycroft Anterior Chamber Cannula, 30 Gauge, Angled 45° 5mm from bend to tip, 22mm tip length



C129

Healon/Viscoelastic


Healon/Viscoelastic Injection Cannula, 21 Gauge, Flattened Round tip, Angled 5mm from bend to tip at 40°, 29mm tip length



C130-25

Rainin

Rainin Air Injection Cannula, 25 Gauge, Angled 7mm from bend to tip, Spatulated tip, 19mm tip length




Cannulas

C130-27

Rainin


Rainin Air Injection Cannula, 27 Gauge, Angled 7mm from bend to tip, Spatulated tip, 19mm tip length



C130-30

Rainin


Rainin Air Injection Cannula, 30 Gauge, Angled 7mm from bend to tip, Spatulated tip, 19mm tip length



C141-5

Gills


Gills Aspirating Cannula, 25 Gauge, Angled 5mm from bend to tip, 18mm tip length



C141-7

Gills


Gills Aspirating Cannula, 25 Gauge, Angled 7mm from bend to tip, 21mm tip length



C141-8

Gills

Gills Aspirating Cannula, 25 Gauge, Angled 8mm from bend to tip, 21mm tip length




Cannulas

C141-9

Gills


Gills Aspirating Cannula, 25 Gauge, Angled 9mm from bend to tip, 21mm tip length



C141-10

Gills


Gills Aspirating Cannula, 25 Gauge, Angled 10mm from bend to tip, 21mm tip length



C143-7

Gills


Gills Aspirating Cannula, 27 Gauge, Angled 7mm from bend to tip, 21mm tip length



C143-10

Gills


Gills Aspirating Cannula, 27 Gauge, Angled 10mm from bend to tip, 21mm tip length



C148

Morrison

Morrison Aspirating Cannula, 23 Gauge, Smooth blunt tip, Angled 7mm from bend to tip at 45°, 25mm tip length




Cannulas

C149

Maumenee


Maumenee Vitreous Aspiring Needle 18 Gauge,
Markings on shaft at 10mm 15mm and 20mm from tip,
25mm tip length



C150

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently
curved shaft, .3mm port, 13mm tip length



C151L

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently
curved shaft with sandblasted tip, .3mm port facing
left, 13mm tip length



C151R

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently
curved shaft with sandblasted tip, .3mm port facing
right, 13mm tip length



C151S

Simcoe

Simcoe Cortex Extractor Cannula, 23 Gauge, Gently
curved shaft with sandblasted tip, .3mm port up
straight, 13mm tip length




Cannulas

C152L

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently curved shaft, .3mm port facing left, 21mm tip length



C152R

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently curved shaft, .3mm port facing right, 21mm tip length



C152S

Simcoe


Simcoe Cortex Extractor Cannula, 23 Gauge, Gently curved shaft, .3mm port up straight, 21mm tip length



C154-21

Simcoe


Simcoe Cortex Extractor Cannula, 21 Gauge, Gently curved shaft, .4mm port, 21mm tip length



C154-23

Simcoe

Simcoe Cortex Extractor Cannula, 23 Gauge Thin Wall, Gently curved shaft, .4mm port, 21mm tip length




Cannulas

C154-23S

Simcoe


Simcoe Cortex Extractor, 23 Gauge Thin Wall, Gently curved shaft with sandblasted tip, .4mm port, 21mm tip length



C158-25

Wiles


Wiles Hydrodissection Cannula, 25 Gauge, J-Shaped Cannula, is gently Angled 135° to facilitate Hydrodissection, 21mm tip length



C158-27

Wiles


Wiles Hydrodissection Cannula, 27 Gauge, J-Shaped Cannula, is gently Angled 135° to facilitate Hydrodissection, 21mm tip length



C159L

Moria


Moria Hydrodissection Cannula, 25 Gauge, J-Shaped 135° tip facing Left with gently vaulted Shaft, 20mm tip length



C159R

Moria

Moria Hydrodissection Cannula, 25 Gauge, J-Shaped 135° tip facing Right with gently vaulted Shaft, 20mm tip length




Cannulas

C160A.3

Anis


Anis Cortex Aspirating Cannula, 23 Gauge Thin Wall, Angled 1.5mm from bend to tip at 35°, .3mm port, 2mm sandblasted tip, 25mm tip length



C160A.4

Anis


Anis Cortex Aspirating Cannula, 23 Gauge Thin Wall, Angled 1.5mm from bend to tip at 35°, .4mm port, 2mm sandblasted tip, 26mm tip length



C160S.3

Anis


Anis Cortex Aspirating Cannula, 23 Gauge Thin Wall, Straight, .3mm port, 2mm sandblasted tip, 25mm tip length



C160S.4

Anis


Anis Cortex Aspirating Cannula, 23 Gauge Thin Wall, Straight, .4mm port, 2mm sandblasted tip, 26mm tip length



C180L

McIntyre-Binkhorst

McIntyre-Binkhorst Cannula, 26 Gauge, Left, Angled 40° at 13mm from bend to tip, with 2.6mm hook shaped tip, 18mm tip length




Cannulas

C180R

McIntyre-Binkhorst


McIntyre-Binkhorst Cannula, 26 Gauge, Right, Angled 40° at 13mm from bend to tip, with 2.6mm hook shaped tip, 18mm tip length



C180S

McIntyre-Binkhorst


McIntyre-Binkhorst Cannula 26 Gauge, Straight, 2.6mm hook shaped tip, 18mm tip length



C181L

Wilbrandt-Binkhorst


Wilbrandt-Binkhorst Cannula, 23 Gauge, Left, Angled 13mm from bend to tip at 30°, with 2.5mm hook shaped tip, 18mm tip length



C181R

Wilbrandt-Binkhorst


Wilbrandt-Binkhorst Cannula, 23 Gauge, Right, Angled 13mm from bend to tip at 30°, with 2.5mm hook shaped tip, 18mm tip length



C181S

Wilbrandt-Binkhorst

Wilbrandt-Binkhorst Cannula, 23 Gauge, Straight, 2.5mm hook shaped tip, 18mm tip length




Cannulas

C190-19

Bishop-Harmon


Bishop-Harmon Cannula, 19 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C190-21

Bishop-Harmon


Bishop-Harmon Cannula, 21 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C190-22

Bishop-Harmon


Bishop-Harmon Cannula, 22 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C190-23

Bishop-Harmon


Bishop-Harmon Cannula, 23 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C190-25

Bishop-Harmon

Bishop-Harmon Cannula, 25 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length




Cannulas

C190-26

Bishop-Harmon


Bishop-Harmon Cannula, 26 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C190-27

Bishop-Harmon


Bishop-Harmon Cannula, 27 Gauge, Angled 35° at 7mm from bend to spatulated tip, 24mm tip length



C191

Bishop-Harmon


Bishop-Harmon Cannula, 19 Gauge, Angled 30° at 5mm from bend to spatulated tip, 25mm tip length



C197

Willis


Willis Irrigating Cannula, 21 Gauge, Angled 35° at 7mm from bend to .5mm spatulated tip with square corners, 21mm tip length



C198

Fry

Fry Irrigating Cannula, 18 Gauge, Flat spatulated tip, Angled 7mm from bend to tip at 35°, 24mm tip length




Cannulas

C199

Thurmond

Thurmond Nucleus Irrigation Cannula, 19 Gauge, Flat spatulated shaft with beveled tip, Angled 7mm from bend to tip at 35°, 24mm tip length



C200P

Adaptor


Adaptor for Silicone Bulb, Nylon



C205

Silicone Bulb

Replacement Silicone Bulb



C210

Bishop-Harmon


Bishop-Harmon Full Set, 1 Bishop-Harmon 19 Gauge, C130-19, 1 Adaptor for Silicone Bulb C142, 1 Silicone Bulb C141



C211

Park

Park Irrigating Cannula, 19 Gauge, Angled 30° at 12mm from bend to flat spatulated tip, 5 .4mm ports, 32mm tip length




Cannulas

C215-23

Knolle


Knolle Anterior Chamber Irrigating Cannula, 23 Gauge,
Angled 2mm from bend to tip, 19mm tip length



C215-25

Knolle


Knolle Anterior Chamber Irrigating Cannula, 25 Gauge,
Angled 2mm from bend to tip, 19mm tip length



C215-27

Knolle


Knolle Anterior Chamber Irrigating Cannula, 27 Gauge,
Angled 2mm from bend to tip, 19mm tip length



C215-30

Knolle


Knolle Anterior Chamber Irrigating Cannula, 30 Gauge,
Angled 2mm from bend to tip, 19mm tip length



C223-23

Troutman

Troutman Olive tip Cannula, 23 Gauge with 1mm
Stainless Olive-Shaped tip, Gently curved shaft, 27mm
tip length




Cannulas

C223-25

Troutman


Troutman Olive tip Cannula, 25 Gauge, with 1mm Stainless Olive-Shaped tip, Gently curved shaft, 27mm tip length



C223A-25

Troutman


Troutman Olive tip Cannula, 25 Gauge, with 1mm Stainless Olive-Shaped tip, Angled 16mm from bend to tip, 27mm tip length



C230

Gills-Welsh


Gills-Welsh Olive tip 25 Gauge, Angled 10mm from bend to tip at 45°, 1mm stainless ball tip, 26mm tip length



C240L-25

Welsh


Welsh Flat Olive tip Cannula, 25 Gauge, Angled 8mm from bend to tip at 35°, Flattened Olive tip angled left, 25mm tip length



C240R-25

Welsh

Welsh Flat Olive tip Cannula, 25 Gauge, Angled 8mm from bend to tip at 35°, Flattened Olive tip angled right, 25mm tip length




Cannulas

C240S-23

Welsh


Welsh Flat Olive tip Cannula, 23 Gauge, Angled 8mm from bend to tip at 35°, Flattened Olive tip straight, 25mm tip length



C240S-25

Welsh


Welsh Flat Olive tip Cannula, 25 Gauge, Angled 8mm from bend to tip at 35°, Flattened Olive tip straight, 25mm tip length



C244-19

Randolph


Randolph Cyclodialysis Cannula, 19 Gauge, Flattened shaft, Angled/Vaulted 12mm from bend to tip, 27mm tip length



C244-21

Randolph


Randolph Cyclodialysis Cannula, 21 Gauge, Flattened shaft, Angled/Vaulted 12mm from bend to tip, 27mm tip length



C244-23

Randolph

Randolph Cyclodialysis Cannula, 23 Gauge, Flattened shaft, Angled/Vaulted 12mm from bend to tip, 27mm tip length




Cannulas

C245

Randolph


Randolph Cyclodialysis Cannula, 19 Gauge, Flattened shaft, Angled/Vaulted 13mm from bend to tip, 27mm tip length



C251

Randolph

Randolph Cyclodialysis Cannula, 19 Gauge, Vary Flattened shaft, Angled/Vaulted 12mm from bend to tip, 27mm tip length



C275-20.5

Simcoe

Simcoe Irrigating/Aspirating Cannula, 20 Gauge, .5mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length



C275-23.3

Simcoe

Simcoe Irrigating/Aspirating Cannula, 21 Gauge, .3mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length, 1mm Offset



C275-23.4

Simcoe

Simcoe Irrigating/Aspirating Cannula Standard, 23 Gauge Thin Wall, .4mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length



Cannulas

C275R-23.3

Simcoe

Simcoe Irrigating/Aspirating Cannula Reverse 23 Gauge Thin Wall, .3mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length



C275R-23.4

Simcoe


Simcoe Irrigating/Aspirating Cannula Reverse 23 Gauge Thin Wall, .4mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length



C276-21.3

Simcoe

Simcoe Irrigating/Aspirating Cannula, 21 Gauge, .3mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length, 1mm Offset



C276-23.3

Simcoe

Simcoe Irrigating/Aspirating Cannula Standard, 23 Gauge Thin Wall, .3mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length, 1mm Offset



C276R-23.3

Simcoe

Simcoe Irrigating/Aspirating Cannula Reverse 23 Gauge Thin Wall, .3mm top port, Curved, Supplied with Silicone Tubing and Adaptor, 15mm tip length, 1mm Offset




Cannulas

C290

Drews


Drews Irrigating Vectis, 23 Gauge, with One Front Irrigating Port at 12 O'clock and Two Back Irrigating Ports at 10 and 2 O'clock, 38mm tip length



C295

Drews-Knolle


Drews-Knolle Reverse Irrigating Vectis, 23 Gauge Thin Wall, with One Front Irrigating Port at 12 O'clock and Two Back Irrigating Ports at 10 and 2 O'clock, 38mm tip length



C300

Knolle-Pearce


Knolle-Pearce Irrigating Vectis, 23 Gauge Thin Wall, with Three Outer Irrigating ports at 12, 10, and 2 O'clock, 38mm tip length



C304

Sheets


Sheets Irrigating Vectis, 21 Gauge loop with Reinforced 14g Straight shaft, Three Outer Irrigating Ports at 10, 12, and 2 O'clock, NO Serrations on top of Vectis, 38mm tip length



C305

Sheets

Sheets Irrigating Vectis, 21 Gauge loop with Reinforced 14g Straight shaft, Three Outer Irrigating Ports at 10, 12, and 2 O'clock, Serrated on top of Vectis, 38mm tip length



Cannulas

C307

Sheets


Sheets Straight Irrigating Vectis with Serrations, 23 Gauge Thin Wall with reinforced 16 Gauge, Straight shaft, Three Outer Irrigating ports at 12, 10, and 2 O'clock, 35mm tip length



C308

Sheets


Sheets Straight Irrigating Vectis with NO Serrations, 23 Gauge Thin Wall with reinforced 16 Gauge, Straight shaft, Three Outer Irrigating ports at 12, 10, and 2 O'clock, 38mm tip length



C310

Sheets

Sheets Irrigation Vectis, 20 Gauge loop with reinforced 14g, Angled 19mm from bend to tip at 20°, Three Outer Irrigating Ports at 12, 9, and 3 O'clock, Serrated on top of Vectus, 38mm tip length



C311

Anis


Anis Irrigation Vectis 23 Gauge, Heart Shaped Loop with .3mm dual front Irrigation ports, 38mm tip length



C315L

Simcoe

Simcoe Nucleus Delivery Loop 25 Gauge, Left, 4mm wide by 13mm long loop with Serrations, dual .2mm Irrigating ports, 13mm tip length




Cannulas

C315R

Simcoe


Simcoe Nucleus Delivery Loop 25 Gauge, Right, 4mm wide by 13mm long loop with Serrations, dual .2mm irrigating ports, 13mm tip length



C325-23

Diamatrix


Lacrimal Cannula, 23 Gauge, Blunt tip, Straight shaft, 25mm tip length



C325-25

Diamatrix


Lacrimal Cannula, 25 Gauge, Blunt tip, Straight shaft, 25mm tip length



C325-27

Diamatrix


Lacrimal Cannula, 27 Gauge, Blunt tip, Straight shaft, 25mm tip length



C327-23

Bailey

Bailey Lacrimal Cannula, 23 Gauge with 19 Gauge Reinforced shaft, 6mm blunt tip, 14mm tip length




Cannulas

C327-25

Bailey


Bailey Lacrimal Cannula, 25 Gauge with 21 Gauge
Reinforced shaft, 7mm blunt tip, 15mm tip length



C327-27

Bailey


Bailey Lacrimal Cannula, 27 Gauge with 21 Gauge
Reinforced shaft, 6mm blunt tip, 15mm tip length



C328

Weil


Weil Lacrimal Cannula, 21 Gauge, Blunt tip with gently
curved shaft, 34mm tip length



C329-19

Diamatrix


Lacrimal Cannula, 19 Gauge, Blunt tip, Side port with
straight shaft, 28mm tip length



C329-20

Diamatrix

Lacrimal Cannula, 20 Gauge, Blunt tip, Side port with
straight shaft, 28mm tip length




Cannulas

C330-23

Diamatrix


Lacrimal Cannula, 23 Gauge with 19 Gauge Reinforced shaft, Straight 12mm tip, 32mm tip length



C330-25

Diamatrix


Lacrimal Cannula, 25 Gauge with 19 Gauge Reinforced shaft, Straight 12mm tip, 32mm tip length



C330-27

Diamatrix


Lacrimal Cannula, 27 Gauge with 19 Gauge Reinforced shaft, Straight 12mm tip, 32mm tip length



C335-23

Diamatrix


Lacrimal Cannula, 23 Gauge with 19 Gauge Reinforced shaft, Curved 10mm tip, 31mm tip length



C335-25

Diamatrix

Lacrimal Cannula, 25 Gauge with 19 Gauge Reinforced shaft, Curved 10mm tip, 31mm tip length




Cannulas

C335-27

Diamatrix


Lacrimal Cannula, 27 Gauge with 19 Gauge Reinforced shaft, Curved 10mm tip, 31mm tip length



C338

Shahinien


Shahinian Lacrimal Cannula, 23 Gauge Straight shaft with blunt tapered tip, .3mm side port, 11mm tip length



C339

Burnstine


Burnstine Lacrimal Cannula, 23 Gauge Straight shaft with tapered tip, .15mm open front, 13mm tip length



C340C

McIntyre


McIntyre Lacrimal Cannula, 23 Gauge, Gently curved shaft, Closed tip, .3mm dual side ports, 46mm tip length



C340S

McIntyre

McIntyre Lacrimal Cannula, 23 Gauge, Straight shaft, Closed tip, .3mm dual side ports, 46mm tip length




Cannulas

C350C

West


West Lacrimal Cannula, 23 Gauge, Gently curved shaft,
Closed tip, .3mm bottom port, 30mm tip length



C350S

West


West Lacrimal Cannula, 23 Gauge, Straight shaft,
Closed tip, .3mm single side port, 30mm tip length



C360-19

Atkinson


Atkinson Retrobulbar Needle 19 Gauge, with Sharp tip,
35mm tip length



C360-20

Atkinson


Atkinson Retrobulbar Needle 20 Gauge, with Sharp tip,
35mm tip length



C360-23

Atkinson

Atkinson Retrobulbar Needle 23 Gauge, with Sharp tip,
35mm tip length




Cannulas

C360-25

Atkinson


Atkinson Retrobulbar Needle 25 Gauge, with Sharp tip, 35mm tip length



C360-30

Atkinson


Atkinson Retrobulbar Needle 30 Gauge, with Sharp tip, 35mm tip length



C364

Stevens

Stevens Sub-Tenon Anesthesia Cannula, 19 Gauge, Gently curved shaft, Flattened spatulated tip, 27mm tip length



C365

Lewicky

Lewicky Chamber Maintainer 20 Gauge, with grooved shaft, Supplied with Silicone Tubing and Adaptor



C366

Lewicky

Lewicky Chamber Maintainer 23 Gauge, Threaded tip, Supplied with Silicone Tubing and Adaptor



Cannulas

C370-23

Heslin

Heslin Infusion Cannula, 23 Gauge, sandblasted
Beveled tip, Supplied with Silicone Tubing and Adaptor



C370-25

Heslin

Heslin Infusion Cannula, 25 Gauge, sandblasted
Beveled tip, Supplied with Silicone Tubing and Adaptor



C372-2.5

Diamatrix

Diamatrix Infusion Cannula, 20 Gauge, polished 2.5mm
Beveled tip, Supplied with 25cm Silicone Tubing and
Adaptor



C372-4.0

Diamatrix

Diamatrix Infusion Cannula, 20 Gauge, polished 4.0mm
Beveled tip, Supplied with 25cm Silicone Tubing and
Adaptor



C372-6.0

Diamatrix

Diamatrix Infusion Cannula, 20 Gauge, polished 6.0mm
Beveled tip, Supplied with 25cm Silicone Tubing and
Adaptor



Cannulas

C376-20

Blumenthal

Blumenthal Anterior Chamber Maintainer, 20 Gauge, with 3mm threaded tip, Supplied with 25cm Silicone Tubing and Adaptor



C376-22

Blumenthal

Blumenthal Anterior Chamber Maintainer, 22 Gauge, with 3mm threaded tip, Supplied with 25cm Silicone Tubing and Adaptor



C376-25

Blumenthal

Blumenthal Anterior Chamber Maintainer, 25 Gauge, with 3mm threaded tip, Supplied with 25cm Silicone Tubing and Adaptor



C377

Diamatrix


Diamatrix Anterior Chamber Maintainer, 25 Gauge, with 3mm smooth shaft, Supplied with 25cm Silicone Tubing and Adaptor



C385

Simcoe

Simcoe Capsule Scrubber 23 Gauge, Curved shaft with sandblasted tip, .3mm top port, 13mm tip length




Cannulas

C390

Simcoe


Simcoe Posterior Capsule Scrubber 23 Gauge, Curved shaft, Reverse Aperture, .3mm bottom port, 13mm tip length



C400-21

Kratz

Kratz Capsule Scrubber 21 Gauge, Angled 1mm with sandblasted tip, 21mm tip length



C400-22

Kratz


Kratz Capsule Scrubber 22 Gauge, Angled 1mm with sandblasted tip, 21mm tip length



C400-23

Kratz


Kratz Capsule Scrubber 23 Gauge, Angled 1mm with sandblasted tip, 21mm tip length



C400-25

Kratz

Kratz Capsule Scrubber 25 Gauge, Angled 1mm with sandblasted tip, 21mm tip length




Cannulas

C400-27

Kratz


Kratz Capsule Scrubber 27 Gauge, Angled 1mm with sandblasted tip, 21mm tip length



C410-18

Jensen


Jensen Capsule Polisher Cannula, 18 Gauge, Curved shaft, Olive-Shaped tip with very rough sandblasted stainless tip, 28mm tip length



C410-23

Jensen


Jensen Capsule Polisher Cannula, 23 Gauge, Curved shaft, Olive-Shaped tip with very rough sandblasted stainless tip, 28mm tip length



C410-25

Jensen


Jensen Capsule Polisher Cannula, 25 Gauge, Curved shaft, Olive-Shaped tip with very rough sandblasted stainless tip, 28mm tip length



C410-27

Jensen

Jensen Capsule Polisher Cannula, 27 Gauge, Curved shaft, Olive-Shaped tip with very rough sandblasted stainless tip, 28mm tip length




Cannulas

C420-23

Jensen

Jensen Posterior Capsule Polisher Cannula, 23 Gauge, Curved shaft, Olive-Shaped tip with very rough posterior sandblasted stainless tip, 28mm tip length



C420-25

Jensen

Jensen Posterior Capsule Polisher Cannula, 25 Gauge, Curved shaft, Olive-Shaped tip with very rough posterior sandblasted stainless tip, 28mm tip length



C420-27

Jensen


Jensen Posterior Capsule Polisher Cannula, 27 Gauge, Curved shaft, Olive-Shaped tip with very rough posterior sandblasted stainless tip, 28mm tip length



C422

Jensen


Jensen Micro Posterior Capsule Polisher 27 Gauge, .8mm Olive-Shaped tip with very rough posterior sandblasted stainless tip, 28mm tip length



C431-27

Drews

Drews Cystotome 27 Gauge, with Cystotome Sharp Cutting Edge, 14mm tip length




Cannulas

C431-30

Drews


Drews Cystotome 30 Gauge, with Cystotome Sharp Cutting Edge, 14 Millimeters tip length



C433-22

Diamatrix


Cystotome Formed 22 Gauge, with Small Cystotome Sharp Cutting Edge, Angled/Vaulted 12mm from bend to tip, 15mm tip length



C433-25

Diamatrix


Cystotome Formed 25 Gauge, with Small Cystotome Sharp Cutting Edge, Angled/Vaulted 12mm from bend to tip, 15mm tip length



C433-27

Diamatrix


Cystotome Formed 27 Gauge, with Small Cystotome Sharp Cutting Edge, Angled/Vaulted 12mm from bend to tip, 15mm tip length



C434

McIntyre

McIntyre Cystotome, 23 Gauge, with Sharp Cutting Point Edge, Angled 10mm from bend to tip, 15mm tip length




Cannulas

C436

Kratz


Kratz Lens Manipulating Cannula, 30 Gauge, .2mm diameter tip, front opening, 5mm from hub to bend, 12mm vaulted tip, 17mm tip length



C443-25

Diamatrix


Hydrodissection Cannula, 25 Gauge, Flat spatulated tip, Angled 11mm from bend to tip, 22mm tip length



C443-27

Diamatrix


Hydrodissection Cannula, 27 Gauge, Flat spatulated tip, Angled 11mm from bend to tip, 22mm tip length



C443-30

Diamatrix


Hydrodissection Cannula, 30 Gauge, Flat spatulated tip, Angled 11mm from bend to tip, 22mm tip length



C444-2

Akahoshi

Akahoshi II Hydrodissection Cannula, 27 Gauge, Beveled tip, 1mm Curved from bend to tip, 19mm tip length




Cannulas

C445-25

Nichamin


Nichamin Hydrodissection Cannula, 25 Gauge, Flat spatulated tip, Angled/Vaulted 12mm from bend to tip, 22mm tip length



C445-27

Nichamin


Nichamin Hydrodissection Cannula, 27 Gauge, Flat spatulated tip, Angled/Vaulted 12mm from bend to tip, 22mm tip length



C447

Kellan


Kellan Hydrodissection Cannula, 25 Gauge, Gently curved shaft with flat spatulated blunt tip, 18mm tip length



C450-23

Diamatrix


Hydrodissection Cannula, 23 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, 22mm tip length



C450-25

Diamatrix

Hydrodissection Cannula, 25 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, 22mm tip length




Cannulas

C450-27

Diamatrix


Hydrodissection Cannula, 27 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, 22mm tip length



C450-30

Diamatrix


Hydrodissection Cannula, 30 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, 22mm tip length



C450L-25

Diamatrix


Hydrodissection Cannula, 25 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, Left, 22mm tip length



C450R-25

Diamatrix


Hydrodissection Cannula, 25 Gauge, Angled 8mm from bend to tip, Flat spatulated tip, Right, 22mm tip length



C459

Osher

Osher Lens Vacuuming Cannula, 27 Gauge, Straight Shaft with Curved tip, Closed tip with .15mm Upward facing port, 25mm tip length




Cannulas

C461-20

Diamatrix


Irrigating Cannula, 20 Gauge, Straight shaft with blunt tip, 25mm tip length



C461-21

Diamatrix


Irrigating Cannula, 21 Gauge, Straight shaft with blunt tip, 25mm tip length



C461-22

Diamatrix


Irrigating Cannula, 22 Gauge, Straight shaft with blunt tip, 25mm tip length



C461-30

Diamatrix


Irrigating Cannula, 30 Gauge, Straight shaft with blunt tip, 25mm tip length



C463-17

Diamatrix

Irrigating Cannula, 17 Gauge, Straight shaft with blunt tip, 38mm tip length




Cannulas

C463-20

Diamatrix


Irrigating Cannula, 20 Gauge, Straight shaft with blunt tip, 38mm tip length



C463-23

Diamatrix


Irrigating Cannula, 23 Gauge, Straight shaft with blunt tip, 38mm tip length



C463-25

Diamatrix

Irrigating Cannula, 25 Gauge, Straight shaft with blunt tip, 38mm tip length



C470

Bracken

Bracken Irrigating Cannula, 19 Gauge, Slightly curved tip with smooth beveled edge, Front open, 28mm tip length



C471

Bracken

Bracken Anterior Chamber Cannula, 19 Gauge, Curved 90° Flattened beveled tip, Polished finish, 21mm tip length




Cannulas

C475-21

Slade


Slade Lasik Cannula, 21 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length



C475-22

Slade


Slade Lasik Cannula, 22 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length



C475-23

Slade


Slade Lasik Cannula, 23 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length



C475-25

Slade


Slade Lasik Cannula, 25 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length



C475-26

Slade

Slade Lasik Cannula, 26 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length




Cannulas

C475-27

Slade


Slade Lasik Cannula, 27 Gauge, Angled/Vaulted 9mm from bend to tip, Flattened spatulated tip, Front open, 16mm tip length



C488-23

Banaji


Banaji Lasik Cannula, 23 Gauge, Bullet shaped tip with 6 .2mm irrigation ports with smooth finish (4 from side to side and 2 from top to bottom), Angled/Vaulted 9mm from bend to tip at 45°, 17mm tip length



C488-25

Banaji


Banaji Lasik Cannula, 25 Gauge, Bullet shaped tip with 6 .2mm irrigation ports with smooth finish (4 from side to side and 2 from top to bottom), Angled/Vaulted 9mm from bend to tip at 45°, 17mm tip length



C488-25S

Banaji


Banaji Lasik Cannula, 25 Gauge, Bullet shaped tip with 6 .2mm irrigation ports with sandblasted finish (4 from side to side and 2 from top to bottom), Angled/Vaulted 9mm from bend to tip at 45°, 17mm tip length



C490

Carter

Carter Lasik Cannula, 25 Gauge, Angled 8mm with Strong Vaulted shaft, Flat spatulated tip, Front open, 22mm tip length




Cannulas

C500-23

Buratto


Buratto Lasik Cannula, 23 Gauge, Angled 10mm with Vaulted shaft, Semi-Closed tip, One front port and two side ports, 25mm tip length



C500-23S

Buratto


Buratto Lasik Cannula, 23 Gauge, with sandblasted tip, Angled/Vaulted 10mm from bend to tip, Semi-Closed tip, One front port and two side ports, 25mm tip length



C500-25

Buratto


Buratto Lasik Cannula, 25 Gauge, Angled 10mm with Vaulted shaft, Semi-Closed tip, One front port and two side ports, 25mm tip length



C500-30

Buratto


Buratto Lasik Cannula, 30 Gauge, Angled/Vaulted 10mm from bend to tip, Semi-Closed tip, One front port and two side ports, 21mm tip length



C520-23

McIntyre

McIntyre Anterior Chamber Cannula, 23 Gauge, Blunt tip, Straight, 18mm tip length




Cannulas

C520-26

McIntyre


McIntyre Anterior Chamber Cannula, 26 Gauge, Blunt tip, Straight, 18mm tip length



C520-30

McIntyre


McIntyre Anterior Chamber Cannula, 30 Gauge, Blunt tip, Straight, 18mm tip length



C521-21

McIntyre


McIntyre Anterior Chamber Cannula, 21 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length



C521-23

McIntyre


McIntyre Anterior Chamber Cannula, 23 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length



C521-25

McIntyre

McIntyre Anterior Chamber Cannula, 25 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length




Cannulas

C521-26

McIntyre


McIntyre Anterior Chamber Cannula, 26 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length



C521-27

McIntyre


McIntyre Anterior Chamber Cannula, 27 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length



C521-30

McIntyre


McIntyre Anterior Chamber Cannula, 30 Gauge, Blunt tip, Angled 12mm from bend to tip, 18mm tip length



C536-25

Scheie


Scheie Anterior Chamber Cannula, 25 Gauge, Straight blunt tip, 20mm tip length



C536-27

Scheie

Scheie Anterior Chamber Cannula, 27 Gauge, Straight blunt tip, 20mm tip length




Cannulas

C537

Seeley


Seeley Hydrodissection Cannula, 26 Gauge, Angled 12mm from bend to tip at 45°, Flattened 2mm Horizontal tip, 12mm tip length



C541L

Diamatrix


Hydrodissection Cannula, 27 Gauge, 3mm U-Shaped tip, Left, 25mm tip length



C541R

Diamatrix


Hydrodissection Cannula, 27 Gauge, 3mm U-Shaped tip, Right, 25mm tip length



C541S

Diamatrix


Hydrodissection Cannula, 27 Gauge, 3mm U-Shaped tip, Straight, 25mm tip length



C550-25.1

Osher

Osher J Cannula, 25 Gauge, 1.5mm wide x 1.5mm tall J-Shaped shaped tip, 26mm tip length




Cannulas

C550-25.2

Osher


Osher J Cannula, 25 Gauge, 2.0mm wide x 2.5mm tall
J-Shaped shaped tip, 25mm tip length



C550-27.1

Osher


Osher J Cannula, 27 Gauge, 1.5mm wide x 1.5mm tall
J-Shaped shaped tip, 26mm tip length



C550-27.2

Osher


Osher J Cannula, 27 Gauge, 2.0mm wide x 2.5mm tall
J-Shaped shaped tip, 25mm tip length



C551-25

Gimbel


Gimbel Hydrodissection Cannula, 25 Gauge, Smooth
blunt J-Shaped shaped tip, 21mm tip length



C551-27

Gimbel

Gimbel Hydrodissection Cannula, 27 Gauge, Smooth
blunt J-Shaped shaped tip, 21mm tip length




Cannulas

C551-30

Gimbel


Gimbel Hydrodissection Cannula, 30 Gauge, Smooth blunt J-Shaped shaped tip, 21mm tip length



C553-23

Fasanella


Fasanella Lacrimal Cannula, 23 Gauge with 19 Gauge Reinforced shaft, 22mm gently curved shaft blunt tip, 42mm tip length



C553-25

Fasanella


Fasanella Lacrimal Cannula, 25 Gauge with 20 Gauge Reinforced shaft, 22mm gently curved shaft blunt tip, 42mm tip length



C555

Gimbel


Gimbel Irrigation Cannula, 26 Gauge, Angled 1mm to side with blunt tip, Angled 7mm bend to bend, Front open, 17mm tip length



C558-27.1.2

Chang

Chang Nucleus Hydrossector Cannula 27ga, 1.25mm flat standard tip angled 90° with beveled opening, 16mm total length




Cannulas

C558-27.1.5

Chang


Chang Nucleus Hydrossector Cannula, 27 Gauge, 1.5mm flat Standard tip angled 90° with beveled opening, 16mm tip length



C558-27.3.0

Chang


Chang Nucleus Hydrossector Cannula, 27 Gauge, 3.0mm flat long tip angled 90° with beveled opening, 16mm tip length



C558-30.1.5

Chang


Chang Nucleus Hydrossector Cannula, 30 Gauge, 1.5mm flat Standard tip angled 90° with beveled opening, 16mm tip length



C560H

Vidaurri


Vidaurri Lasik Cannula, 25 Gauge with 17 Gauge Reinforced shaft, Angled 16mm from bend to tip, 12mm/25 Gauge fork shaped tip with 8 Irrigating Ports, Horizontal, 22mm tip length



C560V

Vidaurri

Vidaurri Lasik Cannula, 25 Gauge with 17 Gauge Reinforced shaft, Angled 16mm from bend to tip, 12mm/25 Gauge fork shaped tip with 8 Irrigating Ports, Vertical, 22mm tip length




Cannulas

C562

Gimbel


Gimbel Forward Irrigating Cannula, 27 Gauge, Blunt tip, Angled/Vaulted 10mm from bend to tip, 25mm tip length



C602

Lane


Lane Lasik Cannula, 25 Gauge, Angled 10mm from bend to tip, Flattened 2mm tip with Posterior sandblast, 3 .17mm ports (2 side ports and 1 port facing up), 21mm tip length



C606

Guell


Guell Lasik Irrigation Cannula, 25 Gauge, Shaft slightly flattened with 6 .25mm side ports and 1 .17mm front port, Angled/Vaulted 12mm from bend to tip, 25mm tip length



C607

Diamatrix


Capsular Dye Cannula, 27 Gauge, Angled 8mm from bend to tip with Vaulted shaft, .17mm Oval Port, 10mm tip length



C612

Gorovoy

Gorovoy Irrigating Descemet's Stripper Cannula, 23 Gauge, Horizontal U-shaped flat spatulated tip, .3mm dual sideports facing right and left, Angled 11mm from bend to tip, 30mm tip length




Cannulas

C615

Yalon


Yalon Hydrodissection Cannula, 25 Gauge, Angled
8mm from bend to tip with irrigation port on the side,
22mm tip length



C616

Yalon

Yalon Hydrodissection Cannula, 25 Gauge, Angled
8mm from bend to tip with irrigation port on the
bottom, 22mm tip length



CA00

Adaptor

Cannula Adaptors Male to Male Luer Adapter
for connecting Aspirating Cannulas to Aspiration
Tubing, Nylon



CA01

Adaptor

Cannula Adaptors Female to Female Luer Lock Cannula
Adaptor, Nylon



CH00

Diamatrix

Infusion Handle Plastic Infusion Handle, Male to
glFemale connector, 3.5 Inches total length, Plastic




Cannulas

CH01

Diamatrix

Infusion Handle Plastic Infusion Handle, Male to Tubing connector, 7.5cm total length, Plastic



CH05

Diamatrix

Infusion Handle Stainless Infusion Handle, Male to Tubing connector, 11.5cm total length, Stainless



CH06

Diamatrix

Infusion Handle Stainless Infusion Handle, Male to Tubing connector, 7.5cm total length, Stainless



CH07

Diamatrix

Infusion Handle Stainless Infusion Handle with finger flow control, Male to Tubing connector, 11.5cm total length, Stainless



CH50

Diamatrix

Infusion Handle Stainless Infusion Handle, Male to Male connector, 10cm total length, Stainless



Cannulas

CH51

Diamatrix

Infusion Handle Stainless Infusion Handle, Male to Female luer lock connector, 10cm total length, Stainless



CH55

Diamatrix

Infusion Handle Stainless Infusion Handle, Male to Male Connector, 7.5cm total length, Stainless



CH60

Diamatrix

Bimanual Irrigation Infusion Handle, Male to Female connector, 4.25 Inches total length, Titanium



CH61

Diamatrix

Bimanual Aspiration Infusion Handle, Male to Male connector, 4.25 Inches total length, Titanium





Fixations Rings




Fixations Rings

TS13FX7

Thornton


Thornton Fixation Ring, 13mm swivel head with 7mm cutout opening in front, No degree markings on top of ring, Flat handle, 110mm total length, Titanium



TS13FX7-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 7mm cutout opening in front, 10° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS13FX9-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 9mm cutout opening in front, 10° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS13FX11-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 11mm cutout opening in front, 10° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS13RX7

Thornton

Thornton Fixation Ring, 13mm swivel head with 7mm cutout opening in front, No degree markings on top of ring, Round handle, 110mm total length, Titanium




Fixations Rings

TS13RX7-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 7mm cutout opening in front, 10° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS13RX7-15

Thornton


Thornton Fixation Ring, 13mm swivel head with 7mm cutout opening in front, 15° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS13RX9-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 9mm cutout opening in front, 10° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS13RX11-10

Thornton


Thornton Fixation Ring, 13mm swivel head with 11mm cutout opening in front, 10° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS16FX9

Thornton

Thornton Fixation Ring, 16mm swivel head with 9mm cutout opening in front, No degree markings on top of ring, Flat handle, 110mm total length, Titanium




Fixations Rings

TS16FX9-10

Thornton


Thornton Fixation Ring, 16mm swivel head with 9mm cutout opening in front, 10° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS16FX9-15

Thornton


Thornton Fixation Ring, 16mm swivel head with 9mm cutout opening in front, 15° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS16FX11-10

Thornton


Thornton Fixation Ring, 16mm swivel head with 11mm cutout opening in front, 10° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS16FX11-15

Thornton


Thornton Fixation Ring, 16mm swivel head with 11mm cutout opening in front, 15° markings visible from top of ring, Flat handle, 110mm total length, Titanium



TS16RX9-10

Thornton

Thornton Fixation Ring, 16mm swivel head with 9mm cutout opening in front, 10° markings visible from top of ring, Round handle, 110mm total length, Titanium




Fixations Rings

TS16RX9-15

Thornton


Thornton Fixation Ring, 16mm swivel head with 9mm cutout opening in front, 15° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS16RX11-10

Thornton


Thornton Fixation Ring, 16mm swivel head with 11mm cutout opening in front, 10° markings visible from top of ring, Round handle, 110mm total length, Titanium



TS16RX11-15

Thornton

Thornton Fixation Ring, 16mm swivel head with 11mm cutout opening in front, 15° markings visible from top of ring, Round handle, 110mm total length, Titanium



T125

Nichamin

Nichamin Fixation Ring, 14mm fixed head with cutout opening in front, Very fine teeth, 10 Degree markings visible from top of ring, Diamond knurled handle, 102mm total length, Stainless





Forceps




Forceps

F001

Desmarres

Desmarres Chalazion Forceps with thumb control lock, 24mm solid oval lower plate with 12mm x 20mm open oval upper plate, 92mm total length, Stainless



F002

Francis

Francis Chalazion Forceps with screw style control lock, 14mm solid oval lower plate with 11mm x 14mm open oval upper plate, 92mm total length, Stainless



F003

Hunt

Hunt Chalazion Forceps with thumb control lock, 15mm solid round lower plate with 12mm open round upper plate, Flat serrated handle, 97mm total length, Stainless



F004

Baird

Baird Chalazion Forceps with thumb control lock, 11mm solid round lower plate with 8mm open round upper plate, Flat serrated handle, 97mm total length, Stainless



F040-1.2

Lester

Lester Fixation Forceps, Fine 1x2 teeth with straight shafts, Flat serrated handle, 96mm total length, Stainless



Forceps

F040-2.3

Lester

Lester Fixation Forceps, Fine 2x3 teeth with Straight shafts, Flat serrated handle, 95mm total length, Stainless



F044-5

Green

Green Fixation Forceps, 5mm T-shaped tips with teeth, Flat serrated handle, 100mm total length, Stainless



F044-10

Green

Green Fixation Forceps, 10mm T-shaped tips with teeth, Flat serrated handle, 100mm total length, Stainless



F044-12

Green

Green Fixation Forceps, 12mm T-shaped tips with teeth, Flat serrated handle, 100mm total length, Stainless



F061

Arita

Arita Meibomian Gland Compressor, Angled 9mm from bend to tip, Smooth heavy polished jaws, Flat handle with hole pattern, 106mm total length, Stainless



Forceps

F095

Scleral

Scleral Plug Cross-Action Forceps, Unique Half Moon Shaped Cups, Angled/Vaulted 9mm from bend to tip, Satin finish with polished tip, 105mm total length, Stainless



F201S

Gills-Colibri

Gills-Colibri Utility Forceps, Very fine pointed tips, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat handle with oval cutout, 75mm total length, Stainless



F201T

Gills-Colibri

Gills-Colibri Utility Forceps, Very fine pointed tips, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat handle with oval cutout, 75mm total length, Titanium



F204

Colibri

Colibri Utility Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, flat handle with oval cutout, 75mm total length, Stainless



F205S

Pierse-Colibri

Pierse-Colibri Corneal Forceps, very delicate 1x1 teeth, Angled 2mm from bend to tip set at 45°, .1mm (.25mm hole), Flat handle with oval cutout, 75mm total length, Stainless



Forceps

F205T

Pierse-Colibri

Pierse-Colibri Corneal Forceps, very delicate 1x1 teeth, Angled 2mm from bend to tip set at 45°, .1mm (.25mm hole), Flat handle with oval cutout, 75mm total length, Titanium



F216S

Colibri

Colibri Utility Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45°, Flat handle with oval cutout, 75mm total length, Stainless



F216T

Colibri


Colibri Utility Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45°, Flat handle with oval cutout, 75mm total length, Titanium



F220S

Troutman-Barraquer

Troutman-Barraquer Colibri Corneal Utility Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45°, 6mm tying platform, Flat serrated handle with ball groove, 73mm total length, Stainless



F220T

Troutman-Barraquer

Troutman-Barraquer Colibri Corneal Utility Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45°, 6mm tying platform, Flat serrated handle with ball groove, 73mm total length, Titanium




Forceps

F229S

Troutman-Barraquer

Troutman-Barraquer Colibri Corneal Utility Forceps, .3mm delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, .12mm, Flat serrated handle with ball groove, 73mm total length, Stainless



F229T

Troutman-Barraquer

Troutman-Barraquer Colibri Corneal Utility Forceps, .3mm delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, .12mm, Flat serrated handle with ball groove, 73mm total length, Titanium



F233-1.2

Bonaccolto

Bonaccolto Utility Forceps, 1.2mm x 16mm long platform with straight shafts, Vertical serrated jaws with cross serrated tips, Flat serrated handle with lateral grooves, 108mm total length, Stainless



F233-1.7

Bonaccolto

Bonaccolto Utility Forceps, 1.7mm x 16mm long platform with straight shafts, Vertical serrated jaws with cross serrated tips, Flat serrated handle with lateral grooves, 109mm total length, Stainless



F235S

Harms-Colibri

Harms-Colibri Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat handle with oval cutout, 75mm total length, Stainless



Forceps

F235T

Harms-Colibri

Harms-Colibri Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat handle with oval cutout, 75mm total length, Titanium



F238-1

Colibri

Colibri Forceps, .4mm standard 1x2 teeth, Angled 2mm from bend to tip set at 45° 6mm tying platform, Flat handle with oval cutout, 75mm total length, Stainless



F238-2

Colibri

Colibri Forceps, .4mm standard 1x2 teeth, Angled 2mm from bend to tip set at 45° without tying platform, Flat handle with oval cutout, 75mm total length, Stainless



F240

Castroviejo-Colibri

Castroviejo-Colibri Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° without tying platform, flat handle with hole pattern, 106mm total length, Stainless



F246S

Colibri

Colibri Corneal Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat serrated handle with lateral grooves, 106mm total length, Stainless



Forceps

F246T

Colibri

Colibri Corneal Forceps, .12mm very delicate 1x2 teeth, Angled 2mm from bend to tip set at 45° with 6mm tying platform, Flat serrated handle with lateral grooves, 106mm total length, Titanium



F300S.3

Castroviejo

Castroviejo Suturing Forceps, .3mm delicate 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Stainless



F300S.5

Castroviejo

Castroviejo Suturing Forceps, .5mm heavy 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Stainless



F300S.12

Castroviejo

Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Stainless



F300T.3

Castroviejo

Castroviejo Suturing Forceps, .3mm delicate 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Titanium



Forceps

F300T.5

Castroviejo

Castroviejo Suturing Forceps, .5mm heavy 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Titanium



F300T.12

Castroviejo

Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5mm tying platform, Straight shafts, flat handle with hole pattern, 106mm total length, Titanium



F301S.3

Castroviejo

Castroviejo Suturing Forceps, .3mm delicate 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Stainless



F301S.5

Castroviejo

Castroviejo Suturing Forceps, .5mm heavy 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Stainless



F301S.12

Castroviejo

Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Stainless



Forceps

F301T.3

Castroviejo

Castroviejo Suturing Forceps, .3mm delicate 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Titanium



F301T.5

Castroviejo

Castroviejo Suturing Forceps, .5mm heavy 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Titanium



F301T.12

Castroviejo

Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5.5mm tying platform, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Titanium



F302S

Sanders-Castroviejo

Sanders-Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5.5mm tying platform, Heavy straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Stainless



F302T

Sanders-Castroviejo

Sanders-Castroviejo Suturing Forceps, .12mm very delicate 1x2 teeth, 5.5mm tying platform, Heavy straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Titanium



Forceps

F304S

Osher


Osher Suturing Forceps, .12mm very delicate 1x2 Elschmig teeth, 5.5mm tying platform, Delicate SET for the finest suture, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Stainless



F304T

Osher

Osher Suturing Forceps, .12mm very delicate 1x2 Elschmig teeth, 5.5mm tying platform, Delicate SET for the finest suture, Straight shafts, Flat serrated handle with lateral grooves, 109mm total length, Titanium



F400

Bonn

Bonn Tissue Forceps, .12mm very delicate 1x2 teeth, Straight shafts, flat handle with hole pattern, 85mm total length, Stainless



F401S

Bonn

Bonn Tissue Forceps, .12mm very delicate 1x2 teeth, Straight shafts, Flat serrated handle, 74mm total length, Stainless



F401T

Bonn

Bonn Tissue Forceps, .12mm very delicate 1x2 teeth, Straight shafts, Flat serrated handle, 74mm total length, Titanium



Forceps

F405S

Bonn

Bonn Tissue Forceps, .12mm very delicate 1x2 teeth, 5mm tying platform, Straight shafts, Flat handle with hole pattern, 102mm total length, Stainless



F405T

Bonn

Bonn Tissue Forceps, .12mm very delicate 1x2 teeth, 5mm tying platform, Straight shafts, Flat handle with hole pattern, 102mm total length, Titanium



F445S.3

Bishop-Harmon

Bishop-Harmon Tissue Forceps, .3mm delicate 1x2 teeth, flat handle with hole pattern, 86mm total length, Stainless



F445S.5

Bishop-Harmon

Bishop-Harmon Tissue Forceps, .5mm heavy 1x2 teeth, flat handle with hole pattern, 86mm total length, Stainless



F445S.8

Bishop-Harmon

Bishop-Harmon Tissue Forceps, .8mm extra heavy 1x2 teeth, flat handle with hole pattern, 86mm total length, Stainless



Forceps

F445T.3

Bishop-Harmon

Bishop-Harmon Tissue Forceps, .3mm delicate 1x2 teeth, flat handle with hole pattern, 86mm total length, Titanium



F445T.5

Bishop-Harmon

Bishop-Harmon Tissue Forceps, .5mm heavy 1x2 teeth, flat handle with hole pattern, 86mm total length, Titanium



F550S

Bishop-Harmon

Bishop-Harmon Dressing Forceps, Criss-Cross Serrated tips, Flat handle with hole pattern, 85mm total length, Stainless



F550T

Bishop-Harmon

Bishop-Harmon Dressing Forceps, Criss-Cross Serrated tips, Flat handle with hole pattern, 85mm total length, Titanium



F552-D

Adson

Adson Dressing Forceps, Delicate serrated jaws, Wide flat serrated handle, 125mm total length, Stainless



Forceps

F552-S

Adson

Dressing Forceps, Standard serrated jaws, Wide flat serrated handle, 125mm total length, Stainless



F558S

Fechtner

Fechtner Conjunctiva Forceps, 1mm very delicate Disc-Shaped tip, 6mm tying platform, Flat handle with hole pattern, 102mm total length, Stainless



F558T

Fechtner


Fechtner Conjunctiva Forceps, 1mm very delicate Disc-Shaped tip, 6mm tying platform, Flat handle with hole pattern, 102mm total length, Titanium



F575S

Putterman


Putterman Style Muscle Forceps, Locking Clamp, 12mm Wide Serrated jaws with 6 Cross-Actions Pins, 100mm total length, Stainless



F575T

Putterman

Putterman Style Muscle Forceps, Locking Clamp, 12mm Wide Serrated jaws with 6 Cross-Actions Pins, 100mm total length, Titanium



Forceps

F600S

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Flat handle with hole pattern, 102mm total length, Stainless



F600T

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Flat handle with hole pattern, 102mm total length, Titanium



F608S-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate with line markings at 2.5 and 5.0mm, Angled 12mm from bend to tip with extremely thin triangular grasping tips, flat handle with hole pattern, 90mm total length, Stainless



F608S-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate with line markings at 2.5 and 5.0mm, Angled/Vaulted 12mm from bend to tip with extremely thin triangular grasping tips, flat handle with hole pattern, 90mm total length, Stainless



F608T-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate with line markings at 2.5 and 5.0mm, Angled 12mm from bend to tip with extremely thin triangular grasping tips, flat handle with hole pattern, 90mm total length, Titanium



Forceps

F608T-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate with line markings at 2.5 and 5.0mm, Angled/Vaulted 12mm from bend to tip with extremely thin triangular grasping tips, flat handle with hole pattern, 90mm total length, Titanium



F610S-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Straight knurled handle, Guide pin stop for tip alignment, 105mm total length, Stainless



F610S-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Straight knurled handle, Guide pin stop for tip alignment, 105mm total length, Stainless



F610T-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Straight knurled handle, Guide pin stop for tip alignment, 105mm total length, Titanium



F610T-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Straight knurled handle, Guide pin stop for tip alignment, 105mm total length, Titanium



Forceps

F612S-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Stainless



F612S-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Stainless



F612T-A

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 11mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Titanium



F612T-V

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Titanium



F615S

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle with flat, Guide pin stop for tip alignment, 105mm total length, Stainless



Forceps

F615T

Utrata

Utrata Capsulorhexis Forceps, Very Delicate, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle with flat, Guide pin stop for tip alignment, 105mm total length, Titanium



F620S

Masket-Utrata

Masket-Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Flat handle with hole pattern, 102mm total length, Stainless



F620T

Masket-Utrata

Masket-Utrata Capsulorhexis Forceps, Very Delicate, Angled/Vaulted Shafts 11mm from bend to tip with extremely thin triangular grasping tips, Flat handle with hole pattern, 102mm total length, Titanium



F628S

Gimbel


Gimbel Capsulorhexis Forceps, jaws set to open only 3mm, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Stainless



F628T

Gimbel

Gimbel Capsulorhexis Forceps, jaws set to open only 3mm, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Titanium



Forceps

F630S-A

Lehner-Utrata

Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Stainless



F630S-V

Lehner-Utrata

Lehner-Utrata Capsulorhexis Forceps, Angled/Vaulted Shafts 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Stainless



F630T-A

Lehner-Utrata

Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Titanium



F630T-V

Lehner-Utrata

Lehner-Utrata Capsulorhexis Forceps, Angled/Vaulted Shafts 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Titanium



F634S

O'Gawa-Utrata

O'Gawa-Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Stainless



Forceps

F634T

O'Gawa-Utrata

O'Gawa-Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, Guide pin stop for tip alignment, 105mm total length, Titanium



F635S

Lehner-Utrata

Lehner-Utrata Capsulorhexis Forceps, Angled/Vaulted Shafts 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, 112mm total length, Stainless



F635T

Lehner-Utrata

Lehner-Utrata Capsulorhexis Forceps, Angled/Vaulted Shafts 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled handle, 112mm total length, Titanium



F637S

Utrata

Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip, extremely thin Sharp triangular grasping tips with cross-action serrated platform, Flat serrated handle, 80mm total length, Stainless



F637T

Utrata

Utrata Capsulorhexis Forceps, Angled 12mm from bend to tip, extremely thin Sharp triangular grasping tips with cross-action serrated platform, Flat serrated handle, 80mm total length, Titanium



Forceps

F649S-1.8

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 1.8mm from bend to tip, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Stainless



F649S-2

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 2mm from bend to tip, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Stainless



F649T-1.8

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 1.8mm from bend to tip, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Titanium



F649T-2

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 2mm from bend to tip, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Flat serrated handle, 80mm total length, Titanium



F651

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 1.8mm from bend to tip, Angled 12mm from bend to tip with extremely thin triangular grasping tips with marks on shaft at 5mm and 6mm, Flat serrated handle, 80mm total length, Stainless



Forceps

F655T

Diamatrix

Diamatrix Capsulorhexis Cystotome Scissor-Action Forceps, Angled/Slightly Vaulted 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled with flat handle, 110mm total length, Titanium



F656T

Koch

Koch Capsulorhexis Cystotome Scissor-Action Forceps, Angled 12mm from bend to tip with extremely thin triangular grasping tips, Diamond knurled with flat handle, 110mm total length, Titanium



F660

Utrata

Utrata MICS Capsulorhexis Forceps, jaws set to open only 1.8mm from bend to tip, Angled/Vaulted 12mm from bend to tip with extremely thin triangular grasping tips SHARP, Diamond knurled handle, 105mm total length, Stainless



F661T

Devgan

Devgan Micro Capsulorhexis Cystotome Scissor-Action Forceps with Ultra-Fine tips, jaws set to open only 2mm from bend to tip, Angled/Slightly Vaulted 12mm from bend to tip with Guide lines at 2 and 4mm, Diamond knurled with flat handle, 110mm total length, Titanium



F663S-A

Utrata-Kershner

Utrata-Kershner Capsulorhexis Cystotome Forceps, Angled 11.5mm from bend to tip with extremely thin Sharp triangular grasping tips, Diamond knurled handle with Flat, Guide pin stop for tip alignment, 105mm total length, Stainless



Forceps

F663S-V

Utrata-Kershner

Utrata-Kershner Capsulorhexis Cystotome Forceps, Angled/Vaulted 11.5mm from bend to tip with extremely thin Sharp triangular grasping tips, Diamond knurled handle with Flat, Guide pin stop for tip alignment, 105mm total length, Stainless



F663T-A

Utrata-Kershner


Utrata-Kershner Capsulorhexis Cystotome Forceps, Angled 11.5mm from bend to tip with extremely thin Sharp triangular grasping tips, Diamond knurled handle with Flat, Guide pin stop for tip alignment, 105mm total length, Titanium



F663T-V

Utrata-Kershner

Utrata-Kershner Capsulorhexis Cystotome Forceps, Angled/Vaulted 11.5mm from bend to tip with extremely thin Sharp triangular grasping tips, Diamond knurled handle with Flat, Guide pin stop for tip alignment, 105mm total length, Titanium



F669S

Giannetti

Giannetti MICS Capsulorhexis Forceps, Very delicate, Angled/Vaulted 11mm from bend to tip with Eye-Shaped gap behind tips, Straight knurled handle, Guide pin stop for tip alignment, 101mm total length, Stainless



F669T

Giannetti

Giannetti MICS Capsulorhexis Forceps, Very delicate, Angled/Vaulted 11mm from bend to tip with Eye-Shaped gap behind tips, Straight knurled handle, Guide pin stop for tip alignment, 101mm total length, Titanium



Forceps

F775S

Haberie-McPherson

Haberie-McPherson Haptic Forceps, Angled 8mm from bend to tip with 7.5mm smooth tying platform with diagonal grooves for grasping IOL haptic, flat handle with hole pattern, 102mm total length, Stainless



F775T

Haberie-McPherson

Haberie-McPherson Haptic Forceps, Angled 8mm from bend to tip with 7.5mm smooth tying platform with diagonal grooves for grasping IOL haptic, flat handle with hole pattern, 102mm total length, Titanium



F782

Diamatrix

Diamatrix IOL Lens Loading Forceps, Smooth flat jaws with unique grasping tips, Flat serrated handle with ball groove for control, 107mm total length, Stainless



F800S

McPherson

McPherson Tying Forceps, Straight shafts, 5mm smooth tying platform, flat handle with hole pattern, 102mm total length, Stainless



F800T

McPherson

McPherson Tying Forceps, Straight shafts, 5mm smooth tying platform, flat handle with hole pattern, 102mm total length, Titanium



Forceps

F801S-S

McPherson

McPherson Tying Forceps, Straight shafts, 4mm smooth tying platform, Flat serrated handle, 86mm total length, Stainless



F801S-V

McPherson

McPherson Tying Forceps, Vaulted shafts, 4mm smooth tying platform, Flat serrated handle, 86mm total length, Stainless



F801T-S

McPherson

McPherson Tying Forceps, Straight shafts, 4mm smooth tying platform, Flat serrated handle, 86mm total length, Titanium



F801T-V

McPherson

McPherson Tying Forceps, Vaulted shafts, 4mm smooth tying platform, Flat serrated handle, 86mm total length, Titanium



F805S

Bechert-McPherson

Bechert-McPherson Tying Forceps, Angled 10.5mm from bend to tip at 45°, 10mm smooth tying platform, Flat serrated handle, 84mm total length, Stainless



Forceps

F805T

Bechert-McPherson

Bechert-McPherson Tying Forceps, Angled 10.5mm from bend to tip at 45°, 10mm smooth tying platform, Flat serrated handle, 84mm total length, Titanium



F808S

Sheets-McPherson

Sheets-McPherson Tying Forceps with 1mm of tips is crisscross serrated, Angled 10.5mm from bend to tip at 45°, 10.5mm tying platform, Flat serrated handle, 82mm total length, Stainless



F808T

Sheets-McPherson

Sheets-Kelman-McPherson Tying Forceps with 1mm of tips is crisscross serrated, Angled 10.5mm from bend to tip at 45°, 10.5mm tying platform, Flat serrated handle, 82mm total length, Titanium



F810S-C

Harms-Tubingen

Harms-Tubingen Tying Forceps, Curved shafts, 7mm smooth tying platform, Flat serrated handle, 86mm total length, Stainless



F810S-S

Harms-Tubingen

Harms-Tubingen Tying Forceps, Straight shafts, 6mm smooth tying platform, Flat serrated handle, 86mm total length, Stainless



Forceps

F810T-C

Harms-Tubingen

Harms-Tubingen Tying Forceps, Curved shafts, 7mm smooth tying platform, Flat serrated handle, 86mm total length, Titanium



F810T-S

Harms-Tubingen

Harms-Tubingen Tying Forceps, Straight shafts, 6mm smooth tying platform, Flat serrated handle, 86mm total length, Titanium



F820S-C

Diamatrix

Tying Forceps, Curved shafts, 7mm tying platform, Flat serrated handle with lateral grooves, 108mm total length, Stainless



F820S-S

Diamatrix

Tying Forceps, Straight shafts, 7mm tying platform, Flat serrated handle with lateral grooves, 108mm total length, Stainless



F820T-C

Diamatrix

Tying Forceps, Curved shafts, 7mm tying platform, Flat serrated handle with lateral grooves, 108mm total length, Titanium



Forceps

F820T-S

Diamatrix

Tying Forceps, Straight shafts, 7mm tying platform, Flat serrated handle with lateral grooves, 108mm total length, Titanium



F842S

McPherson

McPherson Tying Forceps, Angled 5.5mm from bend to tip, 5mm smooth tying platform, Flat handle with hole pattern, 100mm total length, Stainless



F842T

McPherson

McPherson Tying Forceps, Angled 5.5mm from bend to tip, 5mm smooth tying platform, Flat handle with hole pattern, 100mm total length, Titanium



F845S

McPherson

McPherson Tying Forceps, Angled 5mm from bend to tip, 4mm smooth tying platform, Flat serrated handle, 85mm total length, Stainless



F845T

McPherson

McPherson Tying Forceps, Angled 5mm from bend to tip, 4mm smooth tying platform, Flat serrated handle, 85mm total length, Titanium



Forceps

F855S

Kelman-McPherson

Kelman-McPherson Tying Forceps, Angled 8mm from bend to tip, 7.5mm smooth tying platform, flat handle with hole pattern, 102mm total length, Stainless



F855T

Kelman-McPherson

Kelman-McPherson Tying Forceps, Angled 8mm from bend to tip, 7.5mm smooth tying platform, flat handle with hole pattern, 102mm total length, Titanium



F857S

Kelman-McPherson

Kelman-McPherson Tying Forceps, Angled 8mm from bend to tip, 5mm smooth tying platform, Flat serrated handle, 84mm total length, Stainless



F857T

Kelman-McPherson

Kelman-McPherson Tying Forceps, Angled 8mm from bend to tip, 5mm smooth tying platform, Flat serrated handle, 84mm total length, Titanium



F865S-C

Tennant

Tennant Tying Forceps, Extra delicate curved smooth jaws, 6mm tying platform, Straight knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



Forceps

F865S-S

Tennant

Tennant Tying Forceps, Extra delicate straight smooth jaws, 6mm tying platform, Straight knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



F865T-C

Tennant

Tennant Tying Forceps, Extra delicate curved smooth jaws, 6mm tying platform, Straight knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F865T-S

Tennant


Tennant Tying Forceps, Extra delicate straight smooth jaws, 6mm tying platform, Straight knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F869S-C

Girard

Girard Ultra-Fine Tying Forceps, Extra delicate curved smooth jaws, 5mm tying platform Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



F869S-S

Girard

Girard Ultra-Fine Tying Forceps, Extra delicate straight smooth jaws, 5mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



Forceps

F869T-C

Girard

Girard Ultra-Fine Tying Forceps, Extra delicate curved smooth jaws, 5mm tying platform Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F869T-S

Girard

Girard Ultra-Fine Tying Forceps, Extra delicate straight smooth jaws, 5mm tying platform Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F872S

Girard

Girard Tying Forceps, Extra delicate straight smooth jaws, 5mm tying platform Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



F872T

Girard


Girard Tying Forceps, Extra delicate straight smooth jaws, 5mm tying platform Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F875S

Tennant

Tennant Tying Forceps, Extra delicate straight smooth jaws, 6mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Stainless



Forceps

F875T

Tennant


Tennant Tying Forceps, Extra delicate straight smooth jaws, 6mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 110mm total length, Titanium



F880S-5

Girard


Girard Tying Forceps, Extra delicate curved smooth jaws with 5mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 108mm total length, Stainless



F880S-6

Tennant

Tennant Tying Forceps, Extra delicate curved smooth jaws with 6mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 108mm total length, Stainless



F880T-5

Girard

Girard Tying Forceps, Extra delicate curved smooth jaws with 5mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 108mm total length, Titanium



F880T-6

Tennant

Tennant Tying Forceps, Extra delicate curved smooth jaws with 6mm tying platform, Diamond knurled handle, Guide pin stop for tip alignment, 108mm total length, Titanium



Forceps

F884S-C

Jaffe

Jaffe Tying Forceps, Curved shafts, Very delicate 6mm smooth tying platform, Set for 8-11 sutures, flat handle with hole pattern, 102mm total length, Stainless



F884S-S

Jaffe

Jaffe Tying Forceps, Straight shafts, Very delicate 6mm smooth tying platform, Set for 8-11 sutures, flat handle with hole pattern, 102mm total length, Stainless



F884T-C

Jaffe

Jaffe Tying Forceps, Curved shafts, Very delicate 6mm smooth tying platform, Set for 8-11 sutures, flat handle with hole pattern, 102mm total length, Titanium



F884T-S

Jaffe

Jaffe Tying Forceps, Straight shafts, Very delicate 6mm smooth tying platform, Set for 8-11 sutures, flat handle with hole pattern, 102mm total length, Titanium



F885S-C

Harms

Harms Tying Forceps, Curved shafts, Standard 6mm smooth tying platform, Set for 7-9 sutures, flat handle with hole pattern, 102mm total length, Stainless



Forceps

F885S-S

Harms

Harms Tying Forceps, Straight shafts, Standard 6mm smooth tying platform, Set for 7-9 sutures, flat handle with hole pattern, 102mm total length, Stainless



F900-1

Jewelers

Jewelers Forceps, #1 Standard straight shafts, Flat smooth handle, 116mm total length, Stainless



F900-2

Jewelers

Jewelers Forceps, #2 Delicate Blunt/Rounded Straight shafts, Flat smooth handle, 116mm total length, Stainless



F900-3

Jewelers

Jewelers Forceps, #3 Delicate straight shafts, Flat smooth handle, 115mm total length, Stainless



F900-3C

Jewelers

Jewelers Forceps, #3c Delicate straight shafts, Flat smooth handle, 110mm total length, Stainless



Forceps

F900-4

Jewelers

Jewelers Forceps, #4 Very delicate straight shafts, Flat smooth handle, 115mm total length, Stainless



F900-5

Jewelers

Jewelers Forceps, #5 Extra delicate straight shafts, Flat smooth handle, 115mm total length, Stainless



F900-5-45

Jewelers

Jewelers Forceps, #5/45 Extra delicate shafts, Angled 4mm from bend to tip at 45°, Flat smooth handle, 114mm total length, Stainless



F900-5-90

Jewelers

Jewelers Forceps, #5/90 Extra delicate shafts, Angled 4mm from bend to tip at 90°, Flat smooth handle, 113mm total length, Stainless



F900-7

Jewelers

Jewelers Forceps, #7 Delicate curved shafts, Flat smooth handle, 115mm total length, Stainless



Forceps

F942S

Buratto

Buratto Lasik Flap Forceps, .8mm very delicate
Disc-Shaped tip with serrated jaws, Vaulted Shafts, flat
handle with hole pattern, 102mm total length, Stainless



F942T


Buratto

Buratto Lasik Flap Forceps, .8mm very delicate
Disc-Shaped tip with serrated jaws, Vaulted Shafts, flat
handle with hole pattern, 102mm total length, Titanium



Handled Instruments





Handled Instruments

H100-A-D

Sinskey


Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H100-A-S

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H100-A-T

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H100-F-A

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series flat handle, 115mm total length



H101-A-D

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless




Handled Instruments

H101-A-S

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H101-A-T

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H101-F-A

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series flat handle, 115mm total length



H100-F-S

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series flat handle, 120mm total length



H100-S-D

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



Handled Instruments

H100-S-S

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Stainless Steel



H100-S-T

Sinskey

Sinskey Lens Manipulating Hook .25mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Titanium



H101-F-S

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series flat handle, 120mm total length



H101-S-D

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H101-S-S

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Stainless Steel



Handled Instruments

H101-S-T

Sinskey


Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Titanium



H105-D

Sinskey


Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip Reversed, Angled 12mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H105-S

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip Reversed, Angled 12mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H105-T

Sinskey

Sinskey II Lens Manipulating Hook .2mm Diameter blunt tip Reversed, Angled 12mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H130-2A-D

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



Handled Instruments

H130-.2A-S

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H130-.2A-T

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H130-.15A-D

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H130-.15A-S

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H130-.15A-T

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H130-.2S-D

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H130-.2S-S

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Stainless Steel



H130-.2S-T

Fenzl

Fenzl Lens Manipulating Hook .2mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Titanium



H130-.15S-D

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H130-.15S-S

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Stainless Steel



Handled Instruments

H130-15S-T

Fenzl

Fenzl Lens Manipulating Hook .15mm Diameter blunt tip, Straight Shaft, Round knurled Handle, 120mm total length, Titanium



H140-A-D

Kuglen


Kuglen Iris Hook Lens Manipulator Push-Pull Model, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H140-A-S

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H140-A-T

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H140-F-A

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series flat handle, 115mm total length



Handled Instruments

H140-F-S

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model,
Straight shaft, Diamatrix Signature Feather-Series flat
handle, 117mm total length



H140-S-D

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model,
Straight shaft, Diamatrix Signature Feather-Series round
handle, 117mm total length, Titanium/Stainless



H140-S-S

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model,
Straight shaft, Round knurled Handle, 117mm total
length, Stainless Steel



H140-S-T

Kuglen

Kuglen Iris Hook Lens Manipulator Push-Pull Model,
Straight shaft, Round knurled Handle, 117mm total
length, Titanium



H170-A-D

Lester

Lester IOL Manipulator .25mm Diameter Teardrop
shaped tip, Angled 10mm from bend to tip, Diamatrix
Signature Feather-Series round handle, 115mm total
length, Titanium/Stainless



Handled Instruments

H170-A-S

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H170-A-T

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H170-S-D

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight shaft, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H170-S-S

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight shaft, Round knurled Handle, 118mm total length, Stainless Steel



H170-S-T

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight shaft, Round knurled Handle, 118mm total length, Titanium



Handled Instruments

H170-V-D

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Angled/Vaulted 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H170-V-S

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H170-V-T

Lester

Lester IOL Manipulator .25mm Diameter Teardrop shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H171-A-D

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H171-A-S

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



Handled Instruments

H171-A-T

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H171-S-D

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H171-S-S

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Straight Shaft, Round knurled Handle, 118mm total length, Stainless Steel



H171-S-T

Lester

Lester II Lens Manipulator .2mm Diameter Teardrop shaped tip, Straight Shaft, Round knurled Handle, 118mm total length, Titanium



H176-D

Maloney

Maloney IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight/Vaulted 7mm from tip, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



Handled Instruments

H176-S

Maloney

Maloney IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight/Vaulted 7mm from tip, Round knurled Handle, 118mm total length, Stainless Steel



H176-T

Maloney

Maloney IOL Manipulator .25mm Diameter Teardrop shaped tip, Straight/Vaulted 7mm from tip, Round knurled Handle, 118mm total length, Titanium



H180-A-D

Osher

Osher .2mm Y-Hooked tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H180-A-S

Osher

Osher .2mm Y-Hooked tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H180-A-T

Osher

Osher .2mm Y-Hooked tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



Handled Instruments

H180-S-D

Osher

Osher .2mm Y-Hooked tip, Straight Shaft, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H180-S-S

Osher

Osher .2mm Y-Hooked tip, Straight Shaft, Round knurled Handle, 118mm total length, Stainless Steel



H180-S-T

Osher

Osher .2mm Y-Hooked tip, Straight Shaft, Round knurled Handle, 118mm total length, Titanium



H185F-1

Graefe

Graefe Strabismus Hook Size # 1, Flattened from side to side, Angled 8mm from bend to tip, Flat smooth handle, 140mm total length, Stainless



H185F-2

Graefe

Graefe Strabismus Hook Size # 2, Flattened from side to side, Angled 9mm from bend to tip, Flat smooth handle, 140mm total length, Stainless



Handled Instruments

H185F-3

Graefe

Graefe Strabismus Hook Size # 3, Flattened from side to side, Angled 10mm from bend to tip, Flat smooth handle, 140mm total length, Stainless



H188F

Gass

Gass Retinal Detachment Hook, Flattened 13mm hook with 1.5mm oval shaped hole, Flat serrated handle, 143mm total length, Stainless



H210-37-D

Tennant

Tennant Nucleus Manipulator .37mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H210-37-S

Tennant

Tennant Nucleus Manipulator .37mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 118mm total length, Stainless Steel



H210-37-T

Tennant

Tennant Nucleus Manipulator .37mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 118mm total length, Titanium



Handled Instruments

H210-6-D

Tennant

Tennant Nucleus Manipulator .6mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless, "Please specify your desired handle design in the special notes when placing your order"



H210-6-S

Tennant

Tennant Nucleus Manipulator .6mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 118mm total length, Stainless Steel, "Please specify your desired handle design in the special notes when placing your order"



H210-6-T

Tennant

Tennant Nucleus Manipulator .6mm Ball Shaped tip, Angled/Vaulted 10mm from bend to tip, Round knurled Handle, 118mm total length, Titanium, "Please specify your desired handle design in the special notes when placing your order"



H215-A-D

Connor

Connor Wand .5mm Ball Shaped tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H215-A-S

Connor

Connor Wand .5mm Ball Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



Handled Instruments

H215-A-T

Connor

Connor Wand .5mm Ball Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H215-C-D

Connor

Connor Wand .5mm Ball Shaped tip, Curved Shaft, Ambidextrous, Diamatrix Signature Feather-Series round handle, 116mm total length, Titanium/Stainless



H215-C-S

Connor

Connor Wand .5mm Ball Shaped tip, Curved Shaft, Ambidextrous, Round knurled Handle, 116mm total length, Stainless Steel



H215-C-T

Connor

Connor Wand .5mm Ball Shaped tip, Curved Shaft, Ambidextrous, Round knurled Handle, 116mm total length, Titanium



H215-S-D

Connor

Connor Wand .5mm Ball Shaped tip, Straight shaft, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



Handled Instruments

H215-S-S

Connor

Connor Wand .5mm Ball Shaped tip, Straight shaft, Round knurled Handle, 117mm total length, Stainless Steel



H215-S-T

Connor

Connor Wand .5mm Ball Shaped tip, Straight shaft, Round knurled Handle, 117mm total length, Titanium



H220

Siebel

Diamatrix Intra Flap Lifter & Retreatment Spatula, Unique double sided instrument, Intra Flap Lifter, Delicate Slightly Pointed tip, Angled 1.5mm from bend to tip, ~ Retreatment Spatula, Unique Z Shaped tip, 12mm Wide Cylindrical Spatula with Slightly Pointed tip, Diamond knurled handle, 106mm total length



H222-D

Machat

Machat PRK/Lasik Spatula 2mm x 5mm Golfclub Shaped Blade with Semi-Sharpe outer edges, Diamatrix Signature Feather-Series round handle, 119mm total length, Titanium/Stainless

No Image Available

H222-S

Machat

Machat PRK/Lasik Spatula, Stainless Steel

No Image Available

Handled Instruments

H222-T

Machat

Machat PRK/Lasik Spatula, Titanium

No Image Available

H226-D

Steinert-Koch

Steinert-Koch Manipulator .75mm wide flat
Teardrop-Shaped tip with V notch cut on tip, Long
vaulted shaft, Diamatrix Signature Feather-Series
round handle, 117mm total length, Titanium/Stainless



H226-S

Steinert-Koch

Steinert-Koch Manipulator .75mm wide flat
Teardrop-Shaped tip with V notch cut on tip,
Long vaulted shaft, Round knurled Handle,
117mm total length, Stainless Steel



H226-T

Steinert-Koch

Steinert-Koch Manipulator .75mm wide flat
Teardrop-Shaped tip with V notch cut on tip, Long
vaulted shaft, Round knurled Handle, 117mm total
length, Titanium



H230-L-D

Lieberman


Lieberman Microfinger & Connor Wand, Unique
double sided instrument, Lieberman Microfinger Unique
C-Shaped Chopper with blunt tip ~ Connor Wand .5mm
Ball Shaped tip, Both Angled 10mm from bend to tip,
Diamatrix Signature Feather-Series round handle, 142mm
total length, Titanium/Stainless



H230-L-S

Lieberman


Lieberman Microfinger & Connor Wand, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ Connor Wand .5mm Ball Shaped tip, Both Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Stainless Steel



H230-L-T

Lieberman

Lieberman Microfinger & Connor Wand, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ Connor Wand .5mm Ball Shaped tip, Both Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Titanium



H230-R-D

Lieberman

Lieberman Microfinger & Connor Wand, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ Connor Wand .5mm Ball Shaped tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 142mm total length, Titanium/Stainless



H230-R-S

Lieberman

Lieberman Microfinger & Connor Wand, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ Connor Wand .5mm Ball Shaped tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 142mm total length, Stainless Steel



H230-R-T

Lieberman

Lieberman Microfinger & Connor Wand, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ Connor Wand .5mm Ball Shaped tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 142mm total length, Titanium



H232-D

Grayson

Grayson Nucleus Manipulator .8mm Flat Olive-Shaped tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handled, 115mm total length, Titanium/Stainless



H232-S

Grayson

Grayson Nucleus Manipulator .8mm Flat Olive-Shaped tip, Angled 10mm from bend to tip, Round knurled Handled, 115mm total length, Stainless Steel



H232-T

Grayson

Grayson Nucleus Manipulator .8mm Flat Olive-Shaped tip, Angled 10mm from bend to tip, Round knurled Handled, 115mm total length, Titanium



H246-D

Femtosecond

Femtosecond DBL-Ended Spatula, 1.0mm Spatula with Beveled edges and Slightly Pointed tip, Straight / 2.0mm Spatula with Beveled edges and Slightly Pointed tip, Slightly Angled at 8mm from bend to tip, Diamatrix Signature Feather-Series round handled, 115mm total length, Titanium/Stainless Please specify your desired handle design in the special notes when placing your order



H246-S

Femtosecond

Femtosecond DBL-Ended Spatula, 1.0mm Spatula with Beveled edges and Slightly Pointed tip, Straight / 2.0mm Spatula with Beveled edges and Slightly Pointed tip, Slightly Angled at 8mm from bend to tip, Round knurled Handled, 115mm total length, Stainless Steel Please specify your desired handle design in the special notes when placing your order



Handled Instruments

H246-T

Femtosecond

Femtosecond DBL-Ended Spatula, 1.0mm Spatula with Beveled edges and Slightly Pointed tip, Straight / 2.0mm Spatula with Beveled edges and Slightly Pointed tip, Slightly Angled at 8mm from bend to tip, Round knurled Handled, 115mm total length, Titanium Please specify your desired handle design in the special notes when placing your order



H263-L-D

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the right, Intended to be used in the left hand, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H263-L-S

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the right, Intended to be used in the left hand, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H263-L-T

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the right, Intended to be used in the left hand, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H263-R-D

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the left, Intended to be used in the right hand, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



Handled Instruments

H263-R-S

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the left, Intended to be used in the right hand, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Stainless Steel



H263-R-T

Haefliger

Haefliger Phaco Cleaver, Triangular Shaped chopper with inferior cutting edge, Tip is pitched to the left, Intended to be used in the right hand, Angled 10mm from bend to tip, Round knurled Handle, 115mm total length, Titanium



H310F

Jaffe

Jaffe Lens Spatula Flat, .5mm wide flat spatula, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series flat handle, 114mm total length



H310R-D

Jaffe

Jaffe Lens Spatula Flat, .5mm wide flat spatula, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 109mm total length, Titanium/Stainless



H310R-S

Jaffe

Jaffe Lens Spatula Flat, .5mm wide flat spatula, Angled 10mm from bend to tip, Round knurled Handle, 109mm total length, Stainless Steel



Handled Instruments

H312-D

Kimura

Kimura Platinum Spatula, Large tear-drop Spatula tip,
Diamatrix Signature Feather-Series round handle

No Image Available

H312-S

Kimura

Kimura Platinum Spatula, Large tear-drop Spatula tip,
Round knurled Handle



H315-8-D

Koch

Koch Phacoemulsification Spatula, .8mm wide gently
curved shaft with Duckbill-Shaped tip, Diamatrix
Signature Feather-Series round handle, 115mm total
length, Titanium/Stainless



H315-8-S

Koch

Koch Phacoemulsification Spatula, .8mm wide gently
curved shaft with Duckbill-Shaped tip, Round knurled
Handle, 115mm total length, Stainless Steel



H315-8N-D

Koch

Koch Phacoemulsification Spatula, .8mm wide gently
curved shaft with Notches on Duckbill-Shaped tip,
Diamatrix Signature Feather-Series round handle,
115mm total length, Titanium/Stainless




Handled Instruments

H315-8N-S

Koch

Koch Phacoemulsification Spatula, .8mm wide gently curved shaft with Notches on Duckbill-Shaped tip, Round knurled Handle, 115mm total length, Stainless Steel



H315-1-D

Koch

Koch Nucleus Spatula, 1mm wide gently curved shaft with Duckbill-Shaped tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H315-1-S

Koch

Koch Nucleus Spatula, 1mm wide gently curved shaft with Duckbill-Shaped tip, Round knurled Handle, 115mm total length, Stainless Steel



H315-1N-D

Koch

Koch Nucleus Spatula, 1mm wide gently curved shaft with Notches on Duckbill-Shaped tip, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H315-1N-S

Koch

Koch Nucleus Spatula, 1mm wide gently curved shaft with Notches on Duckbill-Shaped tip, Round knurled Handle, 115mm total length, Stainless Steel



Handled Instruments

H320-10-D

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H320-10-S

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 10mm from bend to tip, Round knurled Handle, 118mm total length, Stainless Steel



H320-10-T

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 10mm from bend to tip, Round knurled Handle, 118mm total length, Titanium



H320-7-D

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 7mm from bend to tip, Diamatrix Signature Feather-Series round handle, 118mm total length, Titanium/Stainless



H320-7-S

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 7mm from bend to tip, Round knurled Handle, 118mm total length, Stainless Steel



Handled Instruments

H320-7-T

Koch

Koch Stop and Chop Manipulator, .33mm x 1.5mm long tip design, Angled 7mm from bend to tip, Round knurled Handle, 118mm total length, Titanium



H330-D

Drysdale


Drysdale Nucleus Manipulator, Paddle Shaped tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless



H330-S

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



H330-T

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H331-D

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip with Anterior Capsule Polisher, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless




Handled Instruments

H331-S

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip with Anterior Capsule Polisher, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



H331-T

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip with Anterior Capsule Polisher, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H334-D

Drysdale


Drysdale Nucleus Manipulator, Paddle Shaped tip with Notch on front of tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless



H334-S

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip with Notch on front of tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



H334-T

Drysdale

Drysdale Nucleus Manipulator, Paddle Shaped tip with Notch on front of tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H335-D

Drysdale


Drysdale Nucleus Manipulator, Micro Paddle Shaped tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless



H335-S

Drysdale

Drysdale Nucleus Manipulator, Micro Paddle Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



H335-T

Drysdale

Drysdale Nucleus Manipulator, Micro Paddle Shaped tip, Angled 10mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H345-H-D

Pisacano

Pisacano Nucleus Rotator, Horizontal 1mm flattened ball shape tip, Angled 10.5mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless



H345-H-S

Pisacano

Pisacano Nucleus Rotator, Horizontal 1mm flattened ball shape tip, Angled 10.5mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



Handled Instruments

H345-H-T

Pisacano

Pisacano Nucleus Rotator, Horizontal 1mm flattened ball xshape tip, Angled 10.5mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H345-V-D

Pisacano

Pisacano Nucleus Rotator, Vertical 1mm flattened ball shape tip, Angled 10.5mm from bend to tip, Diamatrix Signature Feather-Series round handle, 112mm total length, Titanium/Stainless



H345-V-S

Pisacano

Pisacano Nucleus Rotator, Vertical 1mm flattened ball shape tip, Angled 10.5mm from bend to tip, Round knurled Handle, 112mm total length, Stainless Steel



H345-V-T

Pisacano

Pisacano Nucleus Rotator, Vertical 1mm flattened ball shape tip, Angled 10.5mm from bend to tip, Round knurled Handle, 112mm total length, Titanium



H348-D

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Angled 7mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



Handled Instruments

H348-S

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H348-T

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Titanium



H350-D

Bechert


Bechert Nucleus Rotator, Vertical Y-Shaped tip, Angled 7mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H350-S

Bechert

Bechert Nucleus Rotator, Vertical Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H350-T

Bechert

Bechert Nucleus Rotator, Vertical Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Titanium



Handled Instruments

H351-D

Bechert

Bechert Nucleus Rotator, Narrow Y-Shaped tip, Angled 7mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H351-S

Bechert


Bechert Nucleus Rotator, Narrow Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H351-T

Bechert

Bechert Nucleus Rotator, Narrow Y-Shaped tip, Angled 7mm from bend to tip, Round knurled Handle, 120mm total length, Titanium



H352-D

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Straight shaft, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H352-S

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Straight shaft, Round knurled Handle, 120mm total length, Stainless Steel



Handled Instruments

H352-T

Bechert

Bechert Nucleus Rotator, Y-Shaped tip, Straight shaft, Round knurled Handle, 120mm total length, Titanium



H353-D

Jaffe-Bechert

Jaffe Bechert Nucleus Rotator, Y-Shaped tip, Angled 9mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



H353-S

Jaffe-Bechert

Jaffe Bechert Nucleus Rotator, Y-Shaped tip, Angled 9mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H353-T

Jaffe-Bechert

Jaffe Bechert Nucleus Rotator, Y-Shaped tip, Angled 9mm from bend to tip, Round knurled Handle, 120mm total length, Titanium



H372-D

Maloney

Maloney PRK Spatula, 2mm x 5mm Semi-Sharp spatulated tip, Angled 5mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless



Handled Instruments

H372-S

Maloney

Maloney PRK Spatula, 2mm x 5mm Semi-Sharp spatulated tip, Angled 5mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H380-.25

Castroviejo

Castroviejo Cyclodialysis, Unique double sided instrument, Spatula, .25mm flat posterior surface with Semi-Sharp tip, Angled/Vaulted 10mm from bend to tip / .25mm flat posterior surface with Semi-Sharp tip, Angled/Vaulted 15mm from bend to tip, Flat handle, 135mm total length



H380-.5

Castroviejo

Castroviejo Cyclodialysis, Unique double sided instrument, Spatula, .5mm flat posterior surface with Semi-Sharp tip, Angled/Vaulted 10mm from bend to tip / .5mm flat posterior surface with Semi-Sharp tip, Angled/Vaulted 15mm from bend to tip, Flat handle, 135mm total length



H392-.25-D

Barraquer

Barraquer Iris Spatula, .25mm Cylindrical smooth tip, Angled/Vaulted 15mm from bend to tip, Diamatrix Signature Feather-Series round handle, 122mm total length, Titanium/Stainless



H392-.25-S

Barraquer

Barraquer Iris Spatula, .25mm Cylindrical smooth tip, Angled/Vaulted 15mm from bend to tip, Round knurled Handle, 122mm total length, Stainless Steel



H392-.5-D

Barraquer

Barraquer Iris Spatula, .5mm Cylindrical smooth tip, Angled/Vaulted 15mm from bend to tip, Diamatrix Signature Feather-Series round handle, 122mm total length, Titanium/Stainless



H392-.5-S

Barraquer

Barraquer Iris Spatula, .5mm Cylindrical smooth tip, Angled/Vaulted 15mm from bend to tip, Round knurled Handle, 122mm total length, Stainless Steel



H392-.5H

Barraquer

Barraquer Iris Spatula, .5mm Cylindrical smooth tip, Angled/Vaulted 15mm from bend to tip, Hexagonal Handle, 100mm total length, Stainless



H395-.25-D

Slade

Slade Refractive Spatula, .25mm Diameter Cylindrical smooth blunt tip, Angled/Vaulted 12mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H395-.25-S

Slade

Slade Refractive Spatula, .25mm Diameter Cylindrical smooth blunt tip, Angled/Vaulted 12mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



Handled Instruments

H395-.5-D

Slade

Slade Refractive Spatula, .5mm Diameter Cylindrical smooth blunt tip, Angled/Vaulted 12mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H395-.5-S

Slade

Slade Refractive Spatula, .5mm Diameter Cylindrical smooth blunt tip, Angled/Vaulted 12mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H405-D

Nichamin

Nichamin Triple Chopper, Unique Multi-Directional cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H405-S

Nichamin

Nichamin Triple Chopper, Unique Multi-Directional cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H405-T

Nichamin

Nichamin Triple Chopper, Unique Multi-Directional cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H408-D

Nagahara

Nagahara Phaco Chopper, Centered flat interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H408-S

Nagahara

Nagahara Phaco Chopper, Centered flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H408-T

Nagahara

Nagahara Phaco Chopper, Centered flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H410-L-D

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H410-L-S

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



Handled Instruments

H410-L-T

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H410-R-D

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H410-R-S

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H410-R-T

Nagahara

Nagahara Diamatrix Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H415-60-D

Nagahara

Nagahara Phaco Chopper, Ambidextrous 60° cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 113mm total length, Titanium/Stainless



Handled Instruments

H415-60-S

Nagahara

Nagahara Phaco Chopper, Ambidextrous 60° cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 113mm total length, Stainless Steel



H415-60-T

Nagahara

Nagahara Phaco Chopper, Ambidextrous 60° cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 113mm total length, Titanium



H415-90-D

Nagahara

Nagahara Phaco Chopper, Ambidextrous 90° cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 113mm total length, Titanium/Stainless



H415-90-S

Nagahara

Nagahara Phaco Chopper, Ambidextrous 90° cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 113mm total length, Stainless Steel



H415-90-T

Nagahara

Nagahara Phaco Chopper, Ambidextrous 90° cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 113mm total length, Titanium



Handled Instruments

H420-L-D

Nagahara


Nagahara Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H420-L-S

Nagahara

Nagahara Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H420-L-T

Nagahara

Nagahara Phaco Chopper, Intended to be used in the left hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H420-R-D

Nagahara

Nagahara Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H420-R-S

Nagahara

Nagahara Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



Handled Instruments

H420-R-T

Nagahara

Nagahara Phaco Chopper, Intended to be used in the right hand, Flat interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H425-D

Fukasaku

Fukasaku Nucleus Chopper, Wedge shaped flat interior cutting edge with blunt tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 109mm total length, Titanium/Stainless



H425-S

Fukasaku

Fukasaku Nucleus Chopper, Wedge shaped flat interior cutting edge with blunt tip, Stainless Steel



H425-T

Fukasaku

Fukasaku Nucleus Chopper, Wedge shaped flat interior cutting edge with blunt tip, Titanium



H427-D

Fukasaku

Fukasaku Small Pupil Snapper Hook Wedge shaped flat inferior cutting edge with blunt tip and Front Notch for engaging Iris, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



Handled Instruments

H427-S

Fukasaku

Fukasaku Small Pupil Snapper Hook , Stainless



H427-T

Fukasaku

Fukasaku Small Pupil Snapper Hook, Titanium



H429-D

Mackool

Mackool Phaco Chopper, Micro ball shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H429-S

Mackool

Mackool Phaco Chopper, Micro ball shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H429-T

Mackool

Mackool Phaco Chopper, Micro ball shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



Handled Instruments

H430-45-D

Rosen

Rosen Phaco Splitter 45°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H430-45-S

Rosen

Rosen Phaco Splitter 45°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H430-45-T

Rosen

Rosen Phaco Splitter 45°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H430-60-D

Rosen

Rosen Phaco Splitter 60°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H430-60-S

Rosen

Rosen Phaco Splitter 60°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H430-60-T

Rosen

Rosen Phaco Splitter 60°, Interior cutting edge with blunt tip, Angled 9mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H435-L-D

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge, Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H435-L-S

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 114mm total length, Stainless Steel



H435-L-T

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 114mm total length, Titanium



H435-R-D

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H435-R-S

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 114mm total length, Stainless Steel



H435-R-T

Nichamin

Nichamin Quick Chopper, Eye-Shaped wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 114mm total length, Titanium



H436-L-D

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H436-L-S

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 114mm total length, Stainless Steel



H436-L-T

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 114mm total length, Titanium



Handled Instruments

H436-R-D

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H436-R-S

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 114mm total length, Stainless Steel



H436-R-T

Nichamin

Nichamin Quick Chopper II, Eye-Shaped wedged cutting edge and Slightly Angled tip, Angled/Vaulted 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 114mm total length, Titanium



H439-AL-D

Siebel

Diamatrix Nucleus Chopper, Angled, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



H439-AL-S

Siebel

Diamatrix Nucleus Chopper, Angled, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Round knurled Handle, 117mm total length, Stainless Steel



Handled Instruments

H439-AL-T

Siebel

Diamatrix Nucleus Chopper, Angled, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Round knurled Handle, 117mm total length, Titanium



H439-AR-D

Siebel

Diamatrix Nucleus Chopper, Angled, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the right hand, Slightly Angled 15mm from bend to tip, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



H439-AR-S

Siebel

Diamatrix Nucleus Chopper, Angled, Stainless Steel



H439-AR-T

Siebel

Diamatrix Nucleus Chopper, Angled, Titanium



H439-SL-D

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



Handled Instruments

H439-SL-S

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Round knurled Handle, 117mm total length, Stainless Steel



H439-SL-T

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the left hand, Round knurled Handle, 117mm total length, Titanium



H439-SR-D

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



H439-SR-S

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the right hand, Round knurled Handle, 117mm total length, Stainless Steel



H439-SR-T

Siebel

Diamatrix Nucleus Chopper, Straight, Curved interior cutting edge with .5mm olive shaped tip, Intended to be used in the right hand, Round knurled Handle, 117mm total length, Titanium



Handled Instruments

H440-D

Siebel

Diamatrix Ambidextrous Chopper, Curved interior cutting edge with .5mm olive shaped tip, Ambidextrous, Slightly Angled 15mm from bend to tip, Diamatrix Signature Feather-Series round handle, 117mm total length, Titanium/Stainless



H440-S

Siebel

Diamatrix Ambidextrous Chopper, Curved interior cutting edge with .5mm olive shaped tip, Ambidextrous, Slightly Angled 15mm from bend to tip, Round knurled Handle, 117mm total length, Stainless Steel



H440-T

Siebel

Diamatrix Ambidextrous Chopper, Curved interior cutting edge with .5mm olive shaped tip, Ambidextrous, Slightly Angled 15mm from bend to tip, Round knurled Handle, 117mm total length, Titanium



H442-L-D

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H442-L-S

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 115mm total length, Stainless Steel




Handled Instruments

H442-L-T

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 115mm total length, Titanium



H442-R-D

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H442-R-S

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 115mm total length, Stainless Steel



H442-R-T

Siebel

Diamatrix Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 115mm total length, Titanium



H445-L-D

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless




Handled Instruments

H445-L-S

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled, Stainless Steel



H445-L-T

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled, Titanium



H445-R-D

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 115mm total length, Titanium/Stainless



H445-R-S

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled, Stainless Steel



H445-R-T

Nagahara

Nagahara-Fine Phaco Chopper, Fine interior cutting edge with blunt tip, Angled, Titanium



Handled Instruments

H455-D

Goldberg

Goldberg Sideport Splitter, .8mm, Flat Teardrop-Shaped blunt tip, Angled 10mm from bend to tip, 113mm total length, Titanium/Stainless



H455-S

Goldberg

Goldberg Sideport Splitter, .8mm, Flat Teardrop-Shaped blunt tip, Angled, Stainless Steel



H455-T

Goldberg

Goldberg Sideport Splitter, .8mm, Flat Teardrop-Shaped blunt tip, Angled, Titanium



H481-1

Lewis

Lewis Lens Loop Small Teardrop-Shaped Loop 4x6mm, Flat Handle, 133mm Overall Length, Stainless



H481-2

Lewis

Lewis Lens Loop Large Teardrop-Shaped Loop 5x7mm, Flat Handle, 134mm Overall Length, Stainless



H485-1-D

Steinert

Steinert I Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted 10mm from bend to tip, 118mm total length, Titanium/Stainless



H485-1-S

Steinert

Steinert I Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted, Stainless Steel



H485-1-T

Steinert

Steinert I Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted, Titanium

No Image Available

H485-2-D

Steinert

Steinert II Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted 10mm from bend to tip, 118mm total length, Titanium/Stainless



H485-2-S

Steinert

Steinert II Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted, Stainless Steel



Handled Instruments

H485-2-T

Steinert

Steinert II Nucleus Claw Chopper, J-Shaped chopper with ambidextrous cutting edge, Angled/Vaulted, Titanium



H550-D

Maloney

Maloney Quick-Chop, Unique double sided instrument, Quick-Chop tip, Angled 9mm from bend to tip ~ Micro Sharp Pick tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 133mm total length, Titanium/Stainless



H550-S

Maloney

Maloney Quick-Chop, Unique double sided instrument, Quick-Chop tip, Angled 9mm from bend to tip ~ Micro Sharp Pick tip, Angled 10mm from bend to tip, Round knurled Handle, 133mm total length, Stainless Steel



H550-T

Maloney

Maloney Quick-Chop, Unique double sided instrument, Quick-Chop tip, Angled 9mm from bend to tip ~ Micro Sharp Pick tip, Angled 10mm from bend to tip, Round knurled Handle, 133mm total length, Titanium



H554-L-D

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 143mm total length, Titanium/Stainless



H554-L-S

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 143mm total length, Stainless Steel



H554-L-T

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 143mm total length, Titanium



H554-R-D

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 143mm total length, Titanium/Stainless



H554-R-S

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 143mm total length, Stainless Steel



H554-R-T

Diamatrix

Diamatrix Nucleus Chopper & Safety Quick Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper Curved interior cutting edge with .5mm Olive Shaped tip, Slightly Angled 15mm from bend to tip ~ Safety Quick Chopper, .5mm Olive-Shaped Blade with 1mm bend, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 143mm total length, Titanium



H558-D

Machat

Machat Lasik Retreatment Spatula, Unique double sided instrument, Cylindrical Spatula with Slightly Pointed tip, Angled 8mm from bend to tip ~ .8mm wide spatula, Flat inferior edge with rounded front tip, Angled 1.2mm from bend to tip, Diamatrix Signature Feather-Series round handle, 120mm total length, Titanium/Stainless

No Image Available

H558-S

Machat

Machat Lasik Retreatment Spatula, Unique double sided instrument, Cylindrical Spatula with Slightly Pointed tip, Angled 8mm from bend to tip ~ .8mm wide spatula, Flat inferior edge with rounded front tip, Angled 1.2mm from bend to tip, Round knurled Handle, 120mm total length, Stainless Steel



H558-T

Machat

Machat Lasik Retreatment Spatula, Unique double sided instrument, Cylindrical Spatula with Slightly Pointed tip, Angled 8mm from bend to tip ~ .8mm wide spatula, Flat inferior edge with rounded front tip, Angled 1.2mm from bend to tip, Round knurled Handle, 120mm total length, Titanium



H559-D

Machat

Machat Lasik Corneal Spatula, Unique double sided instrument, Large Cylindrical Spatula with blunt tip, Angled 10mm from bend to tip ~ .6mm wide spatula, Flat inferior edge with rounded front tip, Angled .7mm from bend to tip, Diamatrix Signature Feather-Series round handle, 119mm total length, Titanium/Stainless

No Image Available

H559-S

Machat

Machat Lasik Corneal Spatula, Unique double sided instrument, Large Cylindrical Spatula with blunt tip, Angled 10mm from bend to tip ~ .6mm wide spatula, Flat inferior edge with rounded front tip, Angled .7mm from bend to tip, Round knurled Handle, 119mm total length, Stainless Steel



H559-T

Machat

Machat Lasik Corneal Spatula, Unique double sided instrument, Large Cylindrical Spatula with blunt tip, Angled 10mm from bend to tip ~ .6mm wide spatula, Flat inferior edge with rounded front tip, Angled .7mm from bend to tip, Round knurled Handle, 119mm total length, Titanium



H560-L-D

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 134mm total length, Titanium/Stainless



H560-L-S

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 134mm total length, Stainless Steel



H560-L-T

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 134mm total length, Titanium



H560-R-D

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 134mm total length, Titanium/Stainless



Handled Instruments

H560-R-S

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 134mm total length, Stainless Steel



H560-R-T

Chang

Chang Microfinger & Quick Chop, Unique double sided instrument, Modified Microfinger, Unique C-Shaped Chopper with blunt tip ~ Quick Chop, 1.3mm Pointed Sharp tip, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 134mm total length, Titanium



H562-L-D

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 142mm total length, Titanium/Stainless



H562-L-S

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 142mm total length, Stainless Steel



H562-L-T

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 142mm total length, Titanium



H562-R-D

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 142mm total length, Titanium/Stainless



H562-R-S

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 142mm total length, Stainless Steel



H562-R-T

Lieberman

Lieberman Microfinger & Spatula, Unique double sided instrument, Lieberman Microfinger Unique C-Shaped Chopper with blunt tip ~ .5mm wide flat spatula, Both Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 142mm total length, Titanium



H568-L-D

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip ~ Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handled, 134mm total length, Titanium/Stainless



H568-L-S

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip ~ Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handled, 134mm total length, Stainless Steel



H568-L-T

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip ~ Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the left hand, Round knurled Handled, 134mm total length, Titanium



H568-R-D

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip, & Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handled, 134mm total length, Titanium/Stainless



H568-R-S

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip, & Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handled, 134mm total length, Stainless Steel



H568-R-T

Nagahara/Nichamin

Nagahara Fine Chopper & Nichamin Quick Chopper, Unique double sided instrument, Nagahara Fine Chopper, Fine interior cutting edge with blunt tip, Angled 10mm from bend to tip, & Nichamin Quick Chopper, Eye-Shaped Wedged cutting edge and Angled tip, Angled 10mm from bend to tip, Intended to be used in the right hand, Round knurled Handled, 134mm total length, Titanium



H573-D

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.0mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 126mm total length, Titanium/Stainless



H573-S

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.0mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Round knurled Handle, 126mm total length, Stainless Steel



H573-T

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.0mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Round knurled Handle, 126mm total length, Titanium

No Image Available

H574-D

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.5mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 126mm total length, Titanium/Stainless



H574-S

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.5mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Round knurled Handle, 126mm total length, Stainless Steel



H574-T

Koch/Nagahara

Koch Spatula and Nagahara Chopper, Unique double sided instrument, Koch Spatula, .8mm Wide Gently curved shaft with Duckbill-Shaped tip, / Nagahara Chopper, 1.5mm - The cutting edge is offset 60° to cut toward Phaco tip, Angled 10mm from bend to tip, Round knurled Handle, 126mm total length, Titanium

No Image Available

H582-L-D

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the left hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 142mm total length, Titanium/Stainless



H582-L-S

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the left hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Stainless Steel



H582-L-T

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the left hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Titanium



H582-R-D

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the right hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 142mm total length, Titanium/Stainless



H582-R-S

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the right hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Stainless Steel



H582-R-T

Chang/Diamatrix

Chang Microfinger Chopper & Diamatrix Vertical Safety Chopper, Unique double sided instrument, Chang Microfinger Chopper, Unique C-Shaped Chopper with blunt tip, Intended to be used in the right hand / Vertical Safety Chopper, .5mm Olive-Shaped Blade with 1mm bend, Both sides Angled 10mm from bend to tip, Round knurled Handle, 142mm total length, Titanium



H584-D

Mackool

Mackool/Cylindrical Spatula, Unique double sided instrument, Mackool Phaco Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, / Cylindrical Spatula, .4mm blade, Angled/Semi-Vaulted 11mm from bend to tip, Diamatrix Signature Feather-Series round handle, 114mm total length, Titanium/Stainless



H584-S

Mackool

Mackool/Cylindrical Spatula, Unique double sided instrument, Mackool Phaco Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, / Cylindrical Spatula, .4mm blade, Angled/Semi-Vaulted 11mm from bend to tip, Round knurled Handle, 114mm total length, Stainless Steel



H584-T

Mackool

Mackool/Cylindrical Spatula, Unique double sided instrument, Mackool Phaco Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, / Cylindrical Spatula, .4mm blade, Angled/Semi-Vaulted 11mm from bend to tip, Round knurled Handle, 114mm total length, Titanium



H588-L-D

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Intended to be used in the left hand, Diamatrix Signature Feather-Series round handle, 138mm total length, Titanium/Stainless



H588-L-S

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 138mm total length, Stainless Steel



H588-L-T

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the left hand, Slightly Angled 15mm from bend to tip, Intended to be used in the left hand, Round knurled Handle, 138mm total length, Titanium



H588-R-D

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the right hand, Slightly Angled 15mm from bend to tip, Intended to be used in the right hand, Diamatrix Signature Feather-Series round handle, 138mm total length, Titanium/Stainless



H588-R-S

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the right hand, Slightly Angled 15mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 138mm total length, Stainless Steel



H588-R-T

Mackool

Diamatrix/Mackool Chopper, Unique double sided instrument, Diamatrix Nucleus Chopper, Micro Ball Shaped blunt tip with interior cutting edge, Angled 10mm from bend to tip, ~ Mackool Phaco Chopper, Curved interior cutting edge with .5mm Olive Shaped tip, Intended to be used in the right hand, Slightly Angled 15mm from bend to tip, Intended to be used in the right hand, Round knurled Handle, 138mm total length, Titanium



H595-D

Lubeck/Akahoshi

Nagahara Lubeck Horizontal Chopper/Akahoshi Nucleus Manipulator, Unique double sided instrument, Nucleus Manipulator, Micro-Ball shaped tip for rotating and manipulating the lens, Angled 10mm from bend to tip ~ Nagahara Phaco Chopper 60° cutting edge, Angled 10mm from bend to tip, Diamatrix Signature Feather-Series round handle, 126mm total length, Titanium/Stainless



H595-S

Lubeck/Akahoshi

Nagahara Lubeck Horizontal Chopper/Akahoshi Nucleus Manipulator, Stainless Steel



H595-T

Lubeck/Akahoshi

Nagahara Lubeck Horizontal Chopper/Akahoshi Nucleus Manipulator, Titanium



H700-0

Quickert

Quickert Lacrimal Intubation Probe Size 0 Malleable Stainless with Pointed tip, 142mm total length



H700-00

Quickert

Quickert Lacrimal Intubation Probe Size 00 Malleable Stainless with Pointed tip, 142mm total length



Handled Instruments

H700-1

Quickert

Quickert Lacrimal Intubation Probe Size 1 Malleable
Stainless with Pointed tip, 142mm total length



H700-2

Quickert

Quickert Lacrimal Intubation Probe Size 2 Malleable
Stainless with Pointed tip, 142mm total length



H700-3

Quickert

Quickert Lacrimal Intubation Probe Size 3 Malleable
Stainless with Pointed tip, 142mm total length



H700-4

Quickert


Quickert Lacrimal Intubation Probe Size 4 Malleable
Stainless with Pointed tip, 142mm total length



H700-SET

Quickert

Quickert Lacrimal Intubation Probe Complete Set 1 of
each H800 through H805 with 10 feet of each Silicone
Tubing (.012" ID and .020" ID)



Handled Instruments

H707-0

Bowman

Bowman Double sided Lacrimal Probe Size # 00-0,
Stainless Malleable with blunt tips



H707-00

Bowman

Bowman Double sided Lacrimal Probe Size # 0000-
000, Stainless Malleable with blunt tips



H707-1.2

Bowman

Bowman Double sided Lacrimal Probe Size # 1-2,
Stainless Malleable with blunt tips



H707-3.4

Bowman

Bowman Double sided Lacrimal Probe Size # 3-4,
Stainless Malleable with blunt tips



H707-5.6

Bowman

Bowman Double sided Lacrimal Probe Size # 5-6,
Stainless Malleable with blunt tips



Handled Instruments

H707-7.8

Bowman


Bowman Double sided Lacrimal Probe Size # 7-8,
Stainless Malleable with blunt tips



H707-SET

Bowman

Bowman Double sided Lacrimal Probe Complete Set,
Stainless Malleable with blunt tips, 1 of each of the
following 0000-000, 00-0, 1-2, 3-4, 5-6, and 7-8



H711

Crawford/Benger

Crawford/Benger Probe Retrieving Hook Flattened
J-Shaped tip, Used for retrieving the Crawford/Benger
Intubation probes, Slightly Angled 10mm from bend to
tip with 65mm Shaft, Diamond knurled handle, 130mm
total length, Stainless



H712

Crawford


Crawford Original Retrieving Hook, Large Flattened
J-Shaped tip, Used for retrieving the Crawford/Benger
Intubation probes, long 75mm Shaft, Diamond knurled
handle, 140mm total length, Stainless



H715

Crawford

Crawford Lacrimal Intubation Set, Two .8mm Diameter
Probes with Olive Shaped tip, Supplied with 10 inches
of Silicone Tubing, Stainless (single Use)



Handled Instruments

H716

Catalan


Catalan Lacrimal Intubation Set, Two .7mm Diameter Probes with blunt tip, Supplied with 30.5cm of Silicone Tubing, Stainless (single Use)

No Image Available

H717

Jackson

Jackson Lacrimal Intubation Set, Two .9mm Diameter Probes with blunt tip, Supplied with 30.5cm of Silicone Tubing, Stainless (single Use)



H718

Diamatrix

Diamatrix Pigtail Probe, Unique double sided instrument, 8mm diameter pigtail curved probes with suture eyes just back of blunt ends, Round handle, 143mm total length, Stainless



H730-0

Wilder

Wilder Lacrimal Dilator Size 0 Infant, .25mm Very fine tapered tip, Diamond knurled handle, 100mm total length, Stainless



H730-1

Wilder

Wilder Lacrimal Dilator Size 1, .25mm Fine Tapered tip, 100mm total length, Stainless



Handled Instruments

H730-2

Wilder

Wilder Lacrimal Dilator Size 2, .4mm Medium Tapered tip, 100mm total length, Stainless



H730-3

Wilder

Wilder Lacrimal Dilator Size 3, .6mm Heavy Tapered tip, 100mm total length, Stainless



H734

Ruedemann

Ruedemann Lacrimal Dilator Short 12mm Long Taper with Needle Type tip, Round Smooth Handle, 70mm total length, Stainless



H735

Castroviejo

Castroviejo Lacrimal Dilator, Unique double sided instrument, Lacrimal Dilator, Size 1 Fine Taper / Size 2 Medium Taper, Diamond knurled handle, 120mm total length, Stainless



H736

Ziegler

Ziegler Lacrimal Dilator, Unique double sided instrument, Gentle Curve on entire Instrument, Diamond knurled handle, 134mm total length, Stainless



Handled Instruments

H803-8

Bunge

Bunge Evisceration Spoon, Small 8mm flat spoon shaped blade, Flat Handle, 136mm total length, Stainless



H803-12

Bunge

Bunge Evisceration Spoon, Large 12mm flat spoon shaped blade, Flat Handle, 136mm total length, Stainless



H805-0

Meyerhoefer

Meyerhoefer Chalazion Curette Size #0, 1.5mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless



H805-1

Meyerhoefer

Meyerhoefer Chalazion Curette Size #1, 1.75mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless



H805-1.5

Meyerhoefer

Meyerhoefer Chalazion Curette Size #1.5, 2.00mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless



Handled Instruments

H805-2

Meyerhoefer

Meyerhoefer Chalazion Curette Size #2, 2.25mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless

No Image Available

H805-3

Meyerhoefer

Meyerhoefer Chalazion Curette Size #3, 2.50mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless



H805-4

Meyerhoefer

Meyerhoefer Chalazion Curette Size #4, 3.50mm Diameter Cup with Semi-Sharp Edges, Flat Handle, 125mm total length, Stainless



H815-L

Shepherd

Shepherd Capsule Polishing Curette, Left, 1.75mm Diameter Ring with Semi-Sharp Top and Bottom Edges, Angled 10mm from bend to tip, Diamond knurled handle, 121mm total length, Stainless



H815-LR

Shepherd

Shepherd Capsule Polishing Curette, Unique double sided instrument, both Right and Left, 1.75mm Diameter low profile Ring with Semi-Sharp Top and Bottom Edges, Slightly Angled 10mm from bend to tip, Diamond knurled handle, 145mm total length, Stainless



Handled Instruments

H815-R

Shepherd

Shepherd Capsule Polishing Curette, Right, 1.75mm Diameter Ring with Semi-Sharp Top and Bottom Edges, Angled 10mm from bend to tip, Diamond knurled handle, 121mm total length, Stainless



H815-S

Shepherd

Shepherd Capsule Polishing Curette, Straight, 1.75mm Diameter Ring with Semi-Sharp Top and Bottom Edges, Angled 10mm from bend to tip, Diamond knurled handle, 121mm total length, Stainless



H850-10

Jameson

Jameson Muscle Hook, 10mm Small ?-Shaped shaft with Paddle Shaped tip, Flat Handle, 130mm total length, Stainless



H850-12

Jameson

Jameson Muscle Hook, 12mm Large ?-Shaped shaft with Paddle Shaped tip, Flat Handle, 130mm total length, Stainless



H853

Stevens

Stevens Tenotomy Hook 6mm Hook-Shaped blunt tip, Flat Handle, 120mm total length, Stainless



Handled Instruments

H854

Green

Green Muscle Hook 10mm Wide ?-Shaped tip with Square blunt tip, Flat Handle, 137mm total length, Stainless



H856

Helveston

Helveston Teaser Hook, Angled 6mm from bend to tip at 90°, Blunt tip, Flat Handle, 115mm total length, Stainless



H857-1

Helveston

Helveston Finder Hook, Small 90° hook with tip slightly angled to a pencil shaped tip, Flat Handle, 117mm total length, Stainless



H857-2

Helveston

Helveston Finder Hook, Large 90° hook with tip slightly angled to a pencil shaped tip, Flat Handle, 117mm total length, Stainless



H871-0

Desmarres

Desmarres Lid Retractor Size #0, Child 11mm Wide Thin Solid Eye-Shaped Blade, Flat Handle, 136mm Overall Length, Stainless



Handled Instruments

H871-1

Desmarres

Desmarres Lid Retractor Size #1, Small 13mm Wide
Thin Solid Eye-Shaped Blade, Flat Handle, 136mm
Overall Length, Stainless



H871-2

Desmarres

Desmarres Lid Retractor Size #2, Medium 15mm Wide
Thin Solid Eye-Shaped Blade, Flat Handle, 136mm
Overall Length, Stainless



H871-3

Desmarres


Desmarres Lid Retractor Size #3, Large 17mm Wide
Thin Solid Eye-Shaped Blade, Flat Handle, 136mm
Overall Length, Stainless



H875-A

Graether

Graether Collar Button Iris retractor and capsule
polisher with roughened posterior surface, Angled
10mm from bend to tip, Diamond knurled handle,
115mm total length, Stainless



H875-S

Graether

Graether Collar Button Iris retractor and capsule
polisher with roughened posterior surface, Straight
shaft, Diamond knurled handle, 118mm total length,
Stainless



Handled Instruments

H880

Schocket

Schocket Scleral Depressor, Unique double sided instrument, T-Shaped Hammer Head tip, Slightly Angled 16mm from bend to tip, ~ Half Ball-Shaped tip, Slightly Angled 16mm from bend to tip, Flat Handle with Pocket clip, 137mm total length, Stainless



H985

Golf

Golf Club Spud with V-Shaped Sharp tip, Flat Handle, 117mm total length, Stainless



H995-10S

Beaver

Beaver Chuck Knife, Diamond knurled handle with Collet for holding Surgical Blades, 10cm total length, Stainless



H995-10S-H

Beaver

Beaver Chuck Knife, Hexagonal handle with Collet for holding Surgical Blades, 10cm total length, Stainless



H995-10T

Beaver

Beaver Chuck Knife, Diamond knurled handle with Collet for holding Surgical Blades, 10cm total length, Titanium



H995-7.5S


Beaver

Beaver Chuck Knife, Diamond knurled handle with Collet for holding Surgical Blades, 7.5cm total length, Stainless



Markers






Markers

M01-3

Thornton


Thornton Low-Profile Pin Marker 3mm,
Diamond knurled handle, 118mm total length



M01-3.5

Thornton


Thornton Low-Profile Pin Marker 3.50mm,
Diamond knurled handle, 118mm total length



M01-3.25

Thornton


Thornton Low-Profile Pin Marker 3.25mm,
Diamond knurled handle, 118mm total length



M01-3.75

Thornton


Thornton Low-Profile Pin Marker 3.75mm,
Diamond knurled handle, 118mm total length



M01-4

Thornton

Thornton Low-Profile Pin Marker 4mm,
Diamond knurled handle, 118mm total length




Markers

M01-4.5

Thornton


Thornton Low-Profile Pin Marker 4.50mm,
Diamond knurled handle, 118mm total
length



M01-4.25

Thornton


Thornton Low-Profile Pin Marker 4.25mm,
Diamond knurled handle, 118mm total
length



M01-4.75

Thornton


Thornton Low-Profile Pin Marker 4.75mm,
Diamond knurled handle, 118mm total
length



M01-5

Thornton


Thornton Low-Profile Pin Marker 5mm,
Diamond knurled handle, 118mm total
length



M01-5.5

Thornton

Thornton Low-Profile Pin Marker 5.50mm,
Diamond knurled handle, 118mm total
length




Markers

M01-5.75

Thornton


Thornton Low-Profile Pin Marker 5.75mm,
Diamond knurled handle, 118mm total
length



M01-6

Thornton


Thornton Low-Profile Pin Marker 6mm,
Diamond knurled handle, 118mm total
length



M01-6.5

Thornton


Thornton Low-Profile Pin Marker 6.50mm,
Diamond knurled handle, 118mm total
length



M01-6.25

Thornton


Thornton Low-Profile Pin Marker 6.25mm,
Diamond knurled handle, 118mm total
length



M01-7

Thornton

Thornton Low-Profile Pin Marker 7mm,
Diamond knurled handle, 118mm total
length




Markers

M01-8

Thornton


Thornton Low-Profile Pin Marker 8mm,
Diamond knurled handle, 118mm total
length



M01-9

Thornton


Thornton Low-Profile Pin Marker 9mm,
Diamond knurled handle, 118mm total
length



M01-10

Thornton


Thornton Low-Profile Pin Marker 10mm,
Diamond knurled handle, 118mm total
length



M01-11

Thornton


Thornton Low-Profile Pin Marker 11mm,
Diamond knurled handle, 118mm total
length



M02-3

Hoffer

Hoffer Optical Marker 3mm with Cross Hairs,
Diamond knurled handle, 118mm total length




Markers

M02-3.25

Hoffer


Hoffer Optical Marker 3.25mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-3.5

Hoffer


Hoffer Optical Marker 3.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-3.75

Hoffer


Hoffer Optical Marker 3.75mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-4

Hoffer


Hoffer Optical Marker 4mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-4.25

Hoffer

Hoffer Optical Marker 4.25mm with Cross Hairs,
Diamond knurled handle, 118mm total length




Markers

M02-4.5

Hoffer


Hoffer Optical Marker 4.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-5

Hoffer


Hoffer Optical Marker 5mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-5.25

Hoffer


Hoffer Optical Marker 5.25mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-5.5

Hoffer


Hoffer Optical Marker 5.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-5.75

Hoffer

Hoffer Optical Marker 5.75mm with Cross Hairs,
Diamond knurled handle, 118mm total length




Markers

M02-6

Hoffer


Hoffer Optical Marker 6mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-6.25

Hoffer


Hoffer Optical Marker 6.25mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-6.5

Hoffer


Hoffer Optical Marker 6.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-7

Hoffer


Hoffer Optical Marker 7mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-7.5

Hoffer

Hoffer Optical Marker 7.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length




Markers

M02-8

Hoffer


Hoffer Optical Marker 8mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-8.5

Hoffer


Hoffer Optical Marker 8.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-8.75

Hoffer


Hoffer Optical Marker 8.75mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-9

Hoffer


Hoffer Optical Marker 9mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-9.5

Hoffer

Hoffer Optical Marker 9.50mm with Cross Hairs,
Diamond knurled handle, 118mm total length




Markers

M02-10

Hoffer


Hoffer Optical Marker 10mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M02-11

Hoffer


Hoffer Optical Marker 11mm with Cross Hairs,
Diamond knurled handle, 118mm total length



M03-3

Bores


Bores Plain Optical Zone Marker 3mm,
Diamond knurled handle, 118mm total length



M03-3.25

Bores


Bores Plain Optical Zone Marker 3.25mm,
Diamond knurled handle, 118mm total
length



M03-3.5

Bores

Bores Plain Optical Zone Marker 3.50mm,
Diamond knurled handle, 118mm total
length




Markers

M03-3.75

Bores


Bores Plain Optical Zone Marker 3.75mm,
Diamond knurled handle, 118mm total
length



M03-4

Bores


Bores Plain Optical Zone Marker 4mm,
Diamond knurled handle, 118mm total
length



M03-4.25

Bores


Bores Plain Optical Zone Marker 4.25mm,
Diamond knurled handle, 118mm total
length



M03-4.5

Bores


Bores Plain Optical Zone Marker 4.50mm,
Diamond knurled handle, 118mm total
length



M03-4.75

Bores

Bores Plain Optical Zone Marker 4.75mm,
Diamond knurled handle, 118mm total
length




Markers

M03-5

Bores


Bores Plain Optical Zone Marker 5mm,
Diamond knurled handle, 118mm total
length



M03-5.5

Bores


Bores Plain Optical Zone Marker 5.50mm,
Diamond knurled handle, 118mm total
length



M03-5.75

Bores


Bores Plain Optical Zone Marker 5.75mm,
Diamond knurled handle, 118mm total
length



M03-6

Bores


Bores Plain Optical Zone Marker 6mm,
Diamond knurled handle, 118mm total
length



M03-6.25

Bores

Bores Plain Optical Zone Marker 6.25mm,
Diamond knurled handle, 118mm total
length




Markers

M03-6.5

Bores


Bores Plain Optical Zone Marker 6.50mm,
Diamond knurled handle, 118mm total
length



M03-7

Bores


Bores Plain Optical Zone Marker 7mm,
Diamond knurled handle, 118mm total
length



M03-7.5

Bores


Bores Plain Optical Zone Marker 7.50mm,
Diamond knurled handle, 118mm total
length



M03-8

Bores


Bores Plain Optical Zone Marker 8mm,
Diamond knurled handle, 118mm total length



M03-8.5

Bores

Bores Plain Optical Zone Marker 8.50mm,
Diamond knurled handle, 118mm total length




Markers

M03-8.75

Bores


Bores Plain Optical Zone Marker 8.75mm,
Diamond knurled handle, 118mm total
length



M03-9

Bores


Bores Plain Optical Zone Marker 9mm,
Diamond knurled handle, 118mm total
length



M03-9.5

Bores


Bores Plain Optical Zone Marker 9.50mm,
Diamond knurled handle, 118mm total
length



M03-10

Bores


Bores Plain Optical Zone Marker 10mm,
Diamond knurled handle, 118mm total
length



M03-11

Bores

Bores Plain Optical Zone Marker 11mm,
Diamond knurled handle, 118mm total
length




Markers

M03F-3

Bores


Bores Plain Optical Zone Marker 3mm,
Flat handle, 118mm total length



M03F-3.25

Bores


Bores Plain Optical Zone Marker 3.25mm,
Flat handle, 118mm total length



M03F-3.5

Bores


Bores Plain Optical Zone Marker 3.50mm,
Flat handle, 118mm total length



M03F-3.75

Bores


Bores Plain Optical Zone Marker 3.75mm,
Flat handle, 118mm total length



M03F-4

Bores

Bores Plain Optical Zone Marker 4mm,
Flat handle, 118mm total length




Markers

M03F-4.25

Bores


Bores Plain Optical Zone Marker 4.25mm,
Flat handle, 118mm total length



M03F-4.5

Bores


Bores Plain Optical Zone Marker 4.50mm,
Flat handle, 118mm total length



M03F-5

Bores


Bores Plain Optical Zone Marker 5mm,
Flat handle, 118mm total length



M03F-6

Bores


Bores Plain Optical Zone Marker 6mm,
Flat handle, 118mm total length



M03F-7

Bores

Bores Plain Optical Zone Marker 7mm,
Flat handle, 118mm total length




Markers

M03F-8

Bores


Bores Plain Optical Zone Marker 8mm,
Flat handle, 118mm total length



M03F-8.5

Bores


Bores Plain Optical Zone Marker 8.50mm,
Flat handle, 118mm total length



M03F-9

Bores

Bores Plain Optical Zone Marker 9mm,
Flat handle, 118mm total length



M05-6T

DSAEK

DSAEK Low-Profile Plain Optical Zone Marker
6mm, Diamond knurled handle, 120mm total
length, Titanium



M05-8.5T

DSAEK

DSAEK Low-Profile Plain Optical Zone Marker
8.50mm, Diamond knurled handle, 120mm
total length, Titanium



Markers

M05-8T

DSAEK

DSAEK Low-Profile Plain Optical Zone Marker
8mm, Diamond knurled handle, 120mm total
length, Titanium



M05-9T

DSAEK


DSAEK Low-Profile Plain Optical Zone Marker
9mm, Diamond knurled handle, 120mm
total length, Titanium



M07-7

Camellin


Camellin Style Lasek 7mm Alcohol Well,
Flat Handle, 102mm total length, Stainless



M07-7.5

Camellin


Camellin Style Lasek 7.5mm Alcohol Well,
Flat Handle, 102mm total length, Stainless



M07-8

Camellin

Camellin Style Lasek 8mm Alcohol Well,
Flat Handle, Stainless




Markers

M07-8.5

Camellin


Camellin Style Lasek 8.5mm Alcohol Well,
Flat Handle, Stainless



M07-9

Camellin


Camellin Style Lasek 9mm Alcohol Well,
Flat Handle, Stainless



M07-9.5

Camellin


Camellin Style Lasek 9.5mm Alcohol Well,
Flat Handle, 102mm total length, Stainless



M08-8

Camellin


Camellin Style Lasek 8mm Trephine,
Flat Handle, 102mm total length,
Stainless



M08-9

Camellin

Camellin Style Lasek 9mm Trephine,
Flat Handle, 102mm total length,
Stainless



Markers

M09-8-8.5

Camellin


Camellin Style Lasek, Unique double sided instrument, 8mm Trephine with 8.5mm Alcohol Well, Flat Handle, Stainless



M09-8-9

Camellin


Camellin Style Lasek, Unique double sided instrument, 8mm Trephine with 9mm Trephine, Flat Handle, 115mm total length, Stainless



M09-8.5-9

Camellin


Camellin Style Lasek, Unique double sided instrument, 8.5mm Trephine with 9mm Alcohol Well, Flat Handle, 115mm total length, Stainless



M09-9-9.5

Camellin


Camellin Style Lasek, Unique double sided instrument, 9mm Trephine with 9.5mm Alcohol Well, Flat Handle, Stainless



M15-8T

DSEK

Corneal DSEK 8mm Marker with Cross-Hairs for providing scoring and stripping of Descemet's membrane during the procedure, Diamond knurled handle, 95mm total length, Titanium




Markers

M15-8.5T

DSEK


Corneal DSEK 8.5mm Marker with Cross-Hairs for providing scoring and stripping of Descemet's membrane during the procedure, Diamond knurled handle, 95mm total length, Titanium



M15-8.5-9T

DSEK


Corneal DSEK Marker, Unique double sided instrument, 8.5mm on one side with 9mm on other, Both sides have Cross-Hairs for providing scoring and stripping of Descemet's membrane during the procedure, Diamond knurled handle, 115mm total length, Titanium



M15-9T

DSEK


Corneal DSEK 9mm Marker with Cross-Hairs for providing scoring and stripping of Descemet's membrane during the procedure, Diamond knurled handle, 95mm total length, Titanium



M15-9.5T

DSEK

Corneal DSEK 9.5mm Marker with Cross-Hairs for providing scoring and stripping of Descemet's membrane during the procedure, Diamond knurled handle, 95mm total length, Titanium



M20-4

Osher-Neumann

Osher-Neumann Corneal Marker 4 Blade, Diamond knurled handle, 118mm total length, Stainless



Markers

M20-6

Osher-Neumann

Osher-Neumann Corneal Marker 6 Blade,
Diamond knurled handle, 118mm total
length, Stainless



M20-8

Osher-Neumann

Osher-Neumann Corneal Marker 8 Blade,
Diamond knurled handle, 118mm total
length, Stainless



M20-12

Osher-Neumann


Osher-Neumann Corneal Marker 12 Blade,
Diamond knurled handle, 118mm total
length, Stainless



M20-16

Osher-Neumann

Osher-Neumann Corneal Marker 16 Blade,
Diamond knurled handle, 118mm total
length, Stainless



M23

Bores

Bores Axis Marker, Dual 4.5mm marking
blades, Will work with any 12mm Mendez
Degree Gauge, Diamond knurled handle,
Stainless



Markers

M25S

Mendez

Mendez Degree Gauge Outside Diameter 18mm, Inside Diameter 12mm, Ruler markings on top of ring, Diamond knurled handle, Stainless



M25T

Mendez


Mendez Degree Gauge Outside Diameter 18mm, Inside Diameter 12mm, Ruler markings on top of ring, Titanium



M28T

Mendez


Mendez Beveled Degree Gauge, Works well with our M58, 5° Ruler markings on top of Unique bevel (labeled at 30° increments), Diamond knurled handle, Titanium



M29T

Mendez


Mendez Beveled Degree Gauge, Works well with our M58, 5° Ruler markings on top of Unique bevel (labeled at 10° and 45° increments), Diamond knurled handle, Titanium



M30S

Mendez

Mendez Degree Gauge with malleable shaft, Outside Diameter 16mm, Inside Diameter 12mm, Ruler markings on top of ring in 5° Increments, Diamond knurled handle, 118mm total length, Stainless




Markers

M30T

Mendez


Mendez Degree Gauge with malleable shaft, Outside Diameter 16mm, Inside Diameter 12mm, Ruler markings on top of ring in 5° Increments, Diamond knurled handle, 118mm total length, Titanium



M33S

Wallace


Wallace Mendez Degree Gauge with malleable shaft, Outside Diameter 16mm, Inside Diameter 12mm, Ruler markings on top of ring in 10° Increments, Diamond knurled handle, 118mm total length, Stainless



M33T

Wallace


Wallace Mendez Degree Gauge with malleable shaft, Outside Diameter 16mm, Inside Diameter 12mm, Ruler markings on top of ring in 10° Increments, Diamond knurled handle, 118mm total length, Titanium



M50S

Diamatrix


Diamatrix Beveled Degree Gauge with four alignment notches, 5° Ruler markings on bevel with handle mounted in the 0° position, 104mm total length, Diamond knurled handle, Stainless



M50T

Diamatrix

Diamatrix Beveled Degree Gauge with four alignment notches, 5° Ruler markings on bevel with handle mounted in the 0° position, 104mm total length, Diamond knurled handle, Titanium




Markers

M51S

Toric


Toric IOL Marker for marking two 2.5mm axis marks, for use with Bevel Degree Gauge, Flat Handle, 100mm total length, Stainless



M51T

Toric


Toric IOL Marker for marking two 2.5mm axis marks, for use with Bevel Degree Gauge, Flat Handle, 100mm total length, Titanium



M52S

Limbal


Limbal Relaxing Incision LRI Marker for marking 30°, 60°, & 90° axis marks, for use with Bevel Degree Gauge, Flat Handle, 95mm total length, Stainless



M52T

Limbal


Limbal Relaxing Incision LRI Marker for marking 30°, 60°, & 90° axis marks, for use with Bevel Degree Gauge, Flat Handle, 95mm total length, Titanium



M53S

Diamatrix

Diamatrix Alignment Gauge Marker with three blunt alignment pins to mark the 3, 6, & 9 O'clock positions, Diamond knurled handle, 100mm total length, Stainless




Markers

M53T

Diamatrix


Diamatrix Alignment Gauge Marker with three blunt alignment pins to mark the 3, 6, & 9 O'clock positions, Diamond knurled handle, 100mm total length, Titanium



M55S

Toric


Toric IOL Reference Marker Blades are positioned at 0°, 90°, 180°, Head is angled at 70°, Diamond knurled handle, 98mm total length, Stainless



M55T

Toric


Toric IOL Reference Marker Blades are positioned at 0°, 90°, 180°, Head is angled at 70°, Diamond knurled handle, 98mm total length, Titanium



M56

Toric

Toric IOL Reference Marker with blades on both sides for marking left and right eye, Blades are positioned at 0°, 90°, 180°, Straight shaft, Diamond knurled handle, 101mm total length, Stainless



M58

Nuijts-Solomon

Nuijts-Solomon Toric IOL Axis Marker with 5mm Optical Zone for use with Mendez Degree Gauge, Diamond knurled handle, 98mm total length, Stainless




Markers

M65T

Diamatrix


Intra-Op Toric Axis Marker with four cross-action 1mm marking blades, Diamond knurled handle mounted in the 90° position, 125mm total length, Titanium



M66T

Diamatrix

Intra-Op Toric Axis Marker II with two 2mm marking blades, Relief cuts on the two 180° positions to allow for easy alignment, Diamond knurled handle mounted in the 90° position, 125mm total length, Titanium



M70S

Castroviejo

Castroviejo Caliper, Straight shaft, 0-20mm measuring scale visible from both sides, Fine threaded control screw for maintaining consistent control, Stainless



M70S-C

Castroviejo

Castroviejo Caliper, Curved shaft, 0-20mm measuring scale visible from both sides, Fine threaded control screw for maintaining consistent control, Stainless



M70T

Castroviejo

Castroviejo Caliper, Straight shaft, 0-20mm measuring scale visible from both sides, Fine threaded control screw for maintaining consistent control, Titanium



Markers

M72

Helveston

Helveston Scleral Ruler single Marking Tooth on tip, 15mm Long Ruler, 1mm Markings with Notches in 5mm increments, Angled/Vaulted 16mm from bend to tip, Flat Handle, 118mm total length, Stainless



M75S

Braunstein

Braunstein Fixed Caliper, Unique double sided instrument, 3mm U-Shaped tip on one side with a 3.5mm U-Shaped tip on the other side, Flat Stock, 70mm total length, Stainless



M75T

Braunstein

Braunstein Fixed Caliper, Unique double sided instrument, 3mm U-Shaped tip on one side with a 3.5mm U-Shaped tip on the other side, Flat Stock, 70mm total length, Titanium



M76S

Braunstein

Braunstein Fixed Caliper, Unique double sided instrument, 3.5mm U-Shaped tip on one side with a 4mm U-Shaped tip on the other side, Flat Stock, 70mm total length, Stainless



M76T

Braunstein

Braunstein Fixed Caliper, Unique double sided instrument, 3.5mm U-Shaped tip on one side with a 4mm U-Shaped tip on the other side, Flat Stock, 70mm total length, Titanium





Needle Holders



Needle Holders

N10-C

Barraquer

Barraquer Needle Holder, Fine Delicate, Curved shaft, 10mm smooth jaws with pointed tips, Without Lock, Diamond knurled handle, 135mm total length, Stainless



N10-CL

Barraquer

Barraquer Needle Holder, Fine Delicate, Curved shaft, 10mm smooth jaws with pointed tips, With Lock, Diamond knurled handle, 135mm total length, Stainless



N10-CLT

Barraquer

Barraquer Needle Holder, Fine Delicate, Curved shaft, 10mm smooth Pointed jaws, With Lock, Diamond knurled handle, 135mm total length, Titanium



N10-CT

Barraquer

Barraquer Needle Holder, Fine Delicate, Curved shaft, 10mm smooth Pointed jaws, Without Lock, Diamond knurled handle, 135mm total length, Titanium



N10-S

Barraquer

Barraquer Needle Holder, Fine Delicate, Straight shaft, 10mm smooth jaws with pointed tips, Without Lock, Diamond knurled handle, 135mm total length, Stainless



Needle Holders

N10-SL

Barraquer

Barraquer Needle Holder, Fine Delicate, Straight shaft, 10mm smooth jaws with pointed tips, With Lock, Diamond knurled handle, 135mm total length, Stainless



N10-SLT

Barraquer

Barraquer Needle Holder, Fine Delicate, Straight shaft, 10mm smooth jaws with pointed tips, With Lock, Diamond knurled handle, 135mm total length, Titanium



N10-ST

Barraquer

Barraquer Needle Holder, Fine Delicate, Straight shaft, 10mm smooth jaws with pointed tips, Without Lock, Diamond knurled handle, 135mm total length, Titanium



N12-C

Barraquer

Barraquer Needle Holder, Delicate, Curved 9mm smooth jaws with cutout, Without Lock, Diamond knurled handle, 135mm total length, Stainless



N12-CL

Barraquer

Barraquer Needle Holder, Delicate, Curved 9mm smooth jaws with cutout, With Lock, Diamond knurled handle, 135mm total length, Stainless



Needle Holders

N12-CLT

Barraquer

Barraquer Needle Holder, Delicate, Curved 9mm smooth jaws with cutout, With Lock, Diamond knurled handle, 135mm total length, Titanium



N12-CT

Barraquer

Barraquer Needle Holder, Delicate, Curved 9mm smooth jaws with cutout, Without Lock, Diamond knurled handle, 135mm total length, Titanium



N12-S

Barraquer

Barraquer Needle Holder, Delicate, Straight shaft 9mm smooth jaws with cutout, Without Lock, Diamond knurled handle, 135mm total length, Stainless



N12-SL

Barraquer

Barraquer Needle Holder, Delicate, Straight shaft 9mm smooth jaws with cutout, With Lock, Diamond knurled handle, 135mm total length, Stainless



N14-C

Barraquer

Barraquer Needle Holder, Standard, Curved 10mm smooth jaws, Without Lock, Diamond knurled handle, 135mm total length, Stainless



Needle Holders

N14-CL

Barraquer

Barraquer Needle Holder, Standard, Curved 10mm smooth jaws, With Lock, Diamond knurled handle, 135mm total length, Stainless



N14-S

Barraquer

Barraquer Needle Holder, Standard, Straight shaft 10mm smooth jaws, Without Lock, Diamond knurled handle, 135mm total length, Stainless



N14-SL

Barraquer


Barraquer Needle Holder, Standard, Straight shaft 10mm smooth jaws, With Lock, Diamond knurled handle, 135mm total length, Stainless



N20-C

Castroviejo

Castroviejo Needle Holder, Delicate, Curved 9mm smooth jaws, Without Lock, Flat serrated handle, 140mm total length, Stainless



N20-CL

Castroviejo

Castroviejo Needle Holder, Delicate, Curved 9mm smooth jaws, With Lock, Flat serrated handle, 140mm total length, Stainless



Needle Holders

N20-S

Castroviejo

Castroviejo Needle Holder, Delicate, Straight shaft
9mm smooth jaws, Without Lock, Flat serrated handle,
140mm total length, Stainless



N20-SL

Castroviejo

Castroviejo Needle Holder, Delicate, Straight shaft
9mm smooth jaws, With Lock, Flat serrated handle,
140mm total length, Stainless



N22-C

Castroviejo

Castroviejo Needle Holder, Standard, Curved 11mm
smooth jaws, Without Lock, Flat serrated handle,
140mm total length, Stainless



N22-CL

Castroviejo

Castroviejo Needle Holder, Standard, Curved 11mm
smooth jaws, With Lock, Flat serrated handle, 140mm
total length, Stainless



N22-S

Castroviejo

Castroviejo Needle Holder, Standard, Straight shaft
11mm smooth jaws, Without Lock, Flat serrated
handle, 140mm total length, Stainless




Needle Holders

N22-SL

Castroviejo

Castroviejo Needle Holder, Standard, Straight shaft
11mm smooth jaws, With Lock, Flat serrated handle,
140mm total length, Stainless



N35-C

Anis

Anis Needle Holder, Extra Delicate, Curved 9mm
smooth jaws, Without Lock, Diamond knurled handle,
127mm total length, Stainless



N35-CL

Anis

Anis Needle Holder, Extra Delicate, Curved 9mm
smooth jaws, With Lock, Diamond knurled handle,
127mm total length, Stainless



N35-S

Anis

Anis Needle Holder, Extra Delicate, Straight shaft 9mm
smooth jaws, Without Lock, Diamond knurled handle,
127mm total length, Stainless



N35-SL

Anis

Anis Needle Holder, Extra Delicate, Straight shaft 9mm
smooth jaws, With Lock, Diamond knurled handle,
127mm total length, Stainless





Scissors



Scissors

S001-A

Vannas

Vannas Capsulotomy Scissors, Sharp Pointed tips, Angled 5mm Blades, Serrated Scissor-Action Handle, 90mm total length, Stainless



S001-C

Vannas

Vannas Capsulotomy Scissors, Sharp Pointed tips, Curved 5mm Blades, Serrated Scissor-Action Handle, 90mm total length, Stainless



S001-S

Vannas

Vannas Capsulotomy Scissors, Sharp Pointed tips, Straight 5mm Blades, Serrated Scissor-Action Handle, 90mm total length, Stainless



S003

Gills-Welsh-Vannas

Gills-Welsh-Vannas Micro Scissors, Sharp Pointed tips, Angled 11mm Blades, 86mm total length, Stainless



S004

McPherson- Vannas

McPherson-Vannas Iris Scissors, Extra-Thin 7mm Curved blade with sharp tips, Serrated Scissor-Action Handle, 86mm total length, Stainless



Scissors

S005-A

Gills-Vannas

Gills-Vannas Capsulotomy Scissors, 7mm Angled blade with sharp pointed tips, Serrated Scissor-Action Handle, 90mm total length, Stainless



S005-C

Gills-Vannas

Gills-Vannas Capsulotomy Scissors, 7mm Curved blade with sharp pointed tips, Serrated Scissor-Action Handle, 90mm total length, Stainless



S005-S

Gills-Vannas

Gills-Vannas Capsulotomy Scissors, 7mm Straight blade with sharp pointed tips, Serrated Scissor-Action Handle, 90mm total length, Stainless



S010-A

Gills-Vannas

Gills-Vannas Scissors, Extra-Thin 7mm angled blades with sharp tips, Three-Hole Scissor-Action Handle, 105mm total length, Stainless



S010-C

Diamatrix-Vannas

Diamatrix-Vannas Scissors, Extra-Thin 7mm Curved blade with sharp pointed tips, Three-Hole Scissor-Action Handle, 105mm total length, Stainless



S010-S

Diamatrix-Vannas

Diamatrix-Vannas Scissors, Extra-Thin 7mm straight blades with sharp tips, Three-Hole Scissor-Action Handle, 105mm total length, Stainless



S012-A

Stern-Gills

Stern-Gills Vannas Scissors, Extra-Thin 10mm angled blades with sharp tips, Three-Hole Scissor-Action Handle, 111mm total length, Stainless



S012-C

Stern-Gills

Stern-Gills Vannas Scissors, Extra-Thin 10mm curved blades with sharp tips, Three-Hole Scissor-Action Handle, 111mm total length, Stainless



S012-S

Stern-Gills

Stern-Gills Vannas Scissors, Extra-Thin 10mm Straight blades with sharp tips, Three-Hole Scissor-Action Handle, 111mm total length, Stainless



S100-C

Diamatrix

Diamatrix Iris Scissors, Curved Blades, Sharp tips, 11mm Blade Length, Serrated Scissor-Action Handle, 102mm total length, Stainless



Scissors

S100-S

Diamatrix

Diamatrix Iris Scissors, Straight Blades, Sharp tips, 11mm Blade Length, Serrated Scissor-Action Handle, 102mm total length, Stainless



S101-C

Diamatrix

Diamatrix Iris Scissor, Standard curved blades with sharp tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



S101-S

Diamatrix

Diamatrix Iris Scissor, Standard straight blades with sharp tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



S105

Barraquer

Barraquer Iris Scissors, 8mm from bend to blunt tips, Serrated Squeeze-Action Handle, 54mm total length, Stainless



S106-C

Diamatrix

Diamatrix Iris Scissors, Curved, 25mm from screw to pointed tips, Ring style handle, 105mm total length, Stainless



Scissors

S106-CS

Diamatrix

Diamatrix Iris Scissors, Lightly Curved, 25mm from screw to pointed tips, Ring style handle, 98mm total length, Stainless



S106-S

Diamatrix

Diamatrix Iris Scissors, Straight, 25mm from screw to pointed tips, Ring style handle, 105mm total length, Stainless



S107-C

Diamatrix

Diamatrix Iris Scissors, Curved, 28mm from screw to pointed tips, Large Ring style handle, 100mm total length, Stainless



S107-S

Diamatrix

Diamatrix Iris Scissors, Straight, 28mm from screw to pointed tips, Large Ring style handle, 100mm total length, Stainless



S200-10

Castroviejo

Castroviejo Universal Corneal Scissors, Curved, Blunt tips, Small Blades, 10mm from pivot to tip, Flat Scissor-Action handle with hole pattern, 112mm total length, Stainless



Scissors

S200-13

Castroviejo

Castroviejo Universal Corneal Scissors, Curved, Blunt tips, Medium Blades, 13mm from pivot to tip, Flat Scissor-Action handle with hole pattern, 120mm total length, Stainless



S200-15

Castroviejo


Castroviejo Universal Corneal Scissors, Curved, Blunt tips, Large Blades, 15mm from pivot to tip, Flat Scissor-Action handle with hole pattern, 120mm total length, Stainless



S204

McGuire

McGuire Corneal Scissors Right, Angled/Curved blades, 21mm from pivot to tip, Ring style handle, 104mm total length, Stainless



S300

Castroviejo

Castroviejo Keratoplasty Scissors, Angled to side, Blunt tips Medium Blades, 16mm from pivot to tip, 108mm total length, Stainless



S350

McPherson-Westcott

McPherson-Westcott Conjunctival Scissors, Curved, Blunt tips Small Blades, 10.5mm from pivot to tip, Flat Scissor-Action handle with hole pattern, 120mm total length, Stainless



Scissors

S401

McPherson-Westcott

McPherson-Westcott Stitch Scissors, Curved, Very Sharp Pointed tips, Small Blades 10.5mm from pivot to tip, Flat Scissor-Action handle with hole pattern, 105mm total length, Stainless



S403

Jaffe

Jaffe Stitch Scissors, Medium curved blades with very sharp tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



S405-B

Westcott

Westcott Stitch Scissors, Curved, Blunt tips 19mm from pivot to tip, Serrated Scissor-Action Handle, 115mm total length, Stainless



S405-S

Westcott

Westcott Stitch Scissors, Curved, Very Sharp Pointed tips, 19mm from pivot to tip, Serrated Scissor-Action Handle, 115mm total length, Stainless



S406-C

Diamatrix

Diamatrix Fine Stitch Scissor, Medium curved blades with sharp tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



Scissors

S406-S

Diamatrix


Diamatrix Fine Stitch Scissor, Medium straight blades with sharp tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



S407-13

Diamatrix

Diamatrix Stitch Scissor, 13mm sharp gently curved blades, Ring style handle, 96mm total length, Stainless



S407-22

Diamatrix


Diamatrix Stitch Scissor, 22mm sharp gently curved blades, Ring style handle, 100mm total length, Stainless



S408

Westcott

Westcott Type Stitch Scissors, 22mm blades from screw to sharp pointed tips, Wide Serrated Scissor-Action Handle, 140mm total length, Stainless



S410-C

Spencer

Spencer Stitch Scissors, Curved blades with bottom blade with suture removal pointed tips, Ring style handle, 89mm total length, Stainless



Scissors

S410-S

Spencer

Spencer Stitch Scissors, Straight blades with bottom blade with suture removal pointed tips, Ring style handle, 89mm total length, Stainless



S502-C

Westcott

Westcott Tenotomy Scissors, Curved, Blunt tips Medium Blades, 19mm from pivot to tip, Serrated Scissor-Action Handle, 117mm total length, Stainless



S502-CL

Westcott

Westcott Tenotomy Scissors, Curved for left handed use, Blunt tips Standard Blades, 21mm from pivot to tip, Serrated Scissor-Action Handle, 117mm total length, Stainless



S502-CR

Westcott

Westcott Tenotomy Scissors, Curved for right handed use, Blunt tips Standard Blades, 21mm from pivot to tip, Serrated Scissor-Action Handle, 117mm total length, Stainless



S507-C

Stevens

Stevens Tenotomy Scissors Light Model, Curved, 19mm from screw to blunt tips, Ring style handle, 98mm total length, Stainless



Scissors

S507-S

Stevens

Stevens Tenotomy Scissors Light Model, Straight, 19mm from screw to blunt tips, Ring style handle, 98mm total length, Stainless



S509

Stevens

Stevens Tenotomy Scissors Heavy Model, Straight, 19mm from screw to blunt tips, Ring style handle, 98mm total length, Stainless



S512

Shepherd-Westcott


Shepard-Westcott Curved Tenotomy Scissors Right, 21mm blades from screw to tip with upper blade serrated, Serrated Scissor-Action Handle, 123mm total length, Stainless



S515-C

Stevens


Stevens Tenotomy Scissors Standard Model, Curved, 30mm from screw to blunt tips, Ring style handle, 112mm total length, Stainless



S515-S

Stevens

Stevens Tenotomy Scissors Standard Model, Straight, 30mm from screw to blunt tips, Ring style handle, 112mm total length, Stainless




Scissors

S516-C

Stevens


Stevens Tenotomy Scissors, 19mm curved blades with blunt tips, Ribbon style handle, 96mm total length, Stainless



S516-S

Stevens

Stevens Tenotomy Scissors, 19mm straight blades with blunt tips, Ribbon style handle, 96mm total length, Stainless



S601-C

Diamatrix

Diamatrix Eye Scissors, Curved, 30mm from screw to blunt tips, Ring style handle, 107mm total length, Stainless



S601-S

Diamatrix

Diamatrix Eye Scissors, Straight, 29mm from screw to blunt tips, Ring style handle, 107mm total length, Stainless



S602-C

Diamatrix

Diamatrix Eye Scissors, Curved, 32mm from screw to sharp tips, Ring style handle, 109mm total length, Stainless



Scissors

S602-S

Diamatrix

Diamatrix Eye Scissors, Straight, 32mm from screw to sharp tips, Ring style handle, 109mm total length, Stainless



S605-C

Diamatrix

Diamatrix Eye Scissors, Curved shafts, Standard pointed tips, Ring style handle, 111mm total length, Stainless



S605-S

Diamatrix

Diamatrix Eye Scissors, Straight shafts, Standard pointed tips, Ring style handle, 111mm total length, Stainless



S608-C

Diamatrix

Diamatrix Eye Scissors, Curved shafts, Standard pointed tips, Large Ring style handle, 111mm total length, Stainless



S608-S

Diamatrix

Diamatrix Eye Scissors, Straight shafts, Standard pointed tips, Large Ring style handle, 111mm total length, Stainless



Scissors

S610-C

Diamatrix


Diamatrix Eye Scissors, Curved shafts, Small pointed tips, Ring style handle, 111mm total length, Stainless



S610-S

Diamatrix

Diamatrix Eye Scissors, Straight shafts, Small pointed tips, Ring style handle, 111mm total length, Stainless



S620-C

Knapp

Knapp Strabismus Scissor, 25mm curved blades with blunt tips, Ring style handle, 114mm total length, Stainless



S620-S

Knapp


Knapp Strabismus Scissor, 25mm straight blades with blunt tips, Ring style handle, 114mm total length, Stainless



S625-C

Strabismus

Strabismus Scissors rounded blades, Curved, 29mm from screw to blunt tips, Ring style handle, 109mm total length, Stainless



Scissors

S625-S

Strabismus


Strabismus Scissors rounded blades, Straight, 29mm from screw to blunt tips, Ring style handle, 109mm total length, Stainless



S741

Diamatrix

Diamatrix Utility Scissors, Slightly Curved, 18mm from screw to blunt tips, Ring style handle, 125mm total length, Stainless



S743

Diamatrix

Diamatrix Utility Scissors, Straight blades with blunt top tips and sharp bottom tip, Ring style handle, 114mm total length, Stainless



S750

Mayo

Mayo Dissecting Scissor, Straight blades, Blunt tips, Ring style handle, 135mm total length, Stainless



S752

Tibolt

Tibolt Punctual Scissors, Standard curved blades with blunt tips, Three-Hole Scissor-Action Handle, 108mm total length, Stainless



Scissors

S800-L

Castroviejo

Castroviejo Corneoscleral Scissors, Medium Blades, 9mm of cutting blade, Left, Serrated Scissor-Action Handle, 105mm total length, Stainless



S800-R

Castroviejo

Castroviejo Corneoscleral Scissors, Medium Blades, 9mm of cutting blade, Right, Serrated Scissor-Action Handle, 105mm total length, Stainless



S805

Metzenbaum

Metzenbaum Dissecting Scissor, Unique slender design with gently curved blades, Blunt tips, Ring style handle, 146mm total length, Stainless



S806

Greenberg

Greenberg Dissecting Scissor, Curved blades, Tungsten-Carbide tips, 32mm from pivot to tip, Ring style handle, 115mm total length, Stainless



S810-1

Enucleation

Enucleation Scissors with heavy rounded blades, Lightly Curved, 40mm from screw to blunt tips, Ring style handle, 141mm total length, Stainless



Scissors

S810-2

Enucleation

Enucleation Scissors with heavy rounded blades, Medium Curved, 44mm from screw to blunt tips, Ring style handle, 139mm total length, Stainless



S810-3

Enucleation

Enucleation Scissors with heavy rounded blades, Strongly Curved, 40mm from screw to blunt tips, Ring style handle, 137mm total length, Stainless





Wire Speculum



Wire Speculum

W00-14

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 14mm



W00-15

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 15mm



W00-16

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 16mm



W00-17

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 17mm



W00-18

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 18mm



Wire Speculum

W00-19

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 19mm



W00-20

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 20mm



W00-21

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 21mm



W00-22

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 22mm



W00-23

Flieringa

Flieringa Scleral Fixation Ring, .020" Stainless wire, Size 23mm




Wire Speculum

W00-SET

Flieringa


Flieringa Scleral Fixation Ring SET, .020" Stainless wire, 14-22mm



W20-13

Lieberman


Lieberman Child lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, Unique triple post speculum, 13mm small blades, Stainless



W20-16

Lieberman


Lieberman Adjustable Unique Triple Post Speculum, This design is intended to lay over the temporal, 16mm Large Blades, Stainless



W20T-16

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, Unique triple post speculum, 16mm Large blades, Titanium



W22-R

Slade

Slade Lasik lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 17mm Open round blades, Stainless




Wire Speculum

W22-S

Slade

Slade Lasik lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 17mm Open square blades with flairs, Stainless



W25-10F

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 10mm Open Square Blades with Flairs, Stainless



W25-14F

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Open Square Blades with Flairs, Stainless



W25N-14F

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the nasal, 14mm Open Square Blades with Flairs, Stainless



W25T-14

Lieberman

Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Open Square Blades, Titanium




Wire Speculum

W25TN-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the bridge of the nose, 14mm Open square blades, Titanium



W28-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Open V-Shaped blades with no flair, Stainless



W28N-14F

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the bridge of the nose, 14mm Open V-Shaped blades with no flairs, Stainless



W28T-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Open V-Shaped blades with no flair, Titanium



W34-10

Lieberman

Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 10mm Solid blades, Stainless




Wire Speculum

W34-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Solid blades, Stainless



W34T-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This design is intended to lay over the temporal, 14mm Solid blades, Titanium



W36-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This unique design is intended to lay either over the temporal or the bridge of the nose, 14mm Solid blades, Stainless



W36T-14

Lieberman


Lieberman lid speculum, Screw style adjustment for solid precise placement, This unique design is intended to lay either over the temporal or the bridge of the nose, 14mm Solid blades, Titanium



W48

Varley

Varley Lasik K-wire lid speculum, .039" Heavy stainless wire, This design is intended to lay over the temporal, 14mm Open square blades with modified bent shafts, 30mm Blade spread, 60mm total length, Stainless




Wire Speculum

W60-10

Barraquer


Barraquer Pediatric lid speculum, .030" Stainless wire, This design is intended to lay over the temporal, 10mm Closed round blades, 15mm Blade spread, 30mm total length, Stainless



W60-14

Barraquer


Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the temporal, 14mm Closed round blades, 18mm Blade spread, 39mm total length, Stainless



W60-14H

Barraquer


Barraquer lid speculum, .039" Heavy Stainless wire, This design is intended to lay over the temporal, 14mm Closed round blades, 20mm Blade spread, 40mm total length, Stainless



W60N-14

Kratz-Barraquer


Kratz-Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the bridge of the nose, 14mm Closed round blades, 18mm Blade spread, 39mm total length, Stainless



W61

Alphonso

Alphonso Newborn lid speculum, .035" Stainless Wire, This design is intended to lay over the temporal, 5mm closed square blade slides inside 8mm closed square blade, Dull finish to reduce glare, 27mm Blade spread, 45mm total length, Stainless




Wire Speculum

W63

Barraquer


Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the bridge of the nose, 14mm Closed round blades, 18mm Blade spread, 39mm total length, Stainless



W65-5

Barraquer


Barraquer Pediatric Speculum, .030" Stainless wire, This design is intended to lay over the temporal, 5mm Closed Square Blades, 15mm Blade Spread, 30mm total length, Stainless



W65-8

Barraquer


Barraquer Pediatric Speculum, .030" Stainless wire, This design is intended to lay over the temporal, 8mm Closed Square Blades, 15mm Blade Spread, 30mm total length, Stainless



W65-9

Barraquer


Barraquer Pediatric Speculum, .030" Stainless Wire, This design is intended to lay over the temporal, 9mm Closed Square Blades, 15mm Blade Spread, 30mm total length, Stainless



W65-11

Barraquer

Barraquer Pediatric Speculum, .030" Stainless Wire, This design is intended to lay over the temporal, 11mm Closed Square Blades, 19mm Blade Spread, 34mm total length, Stainless




Wire Speculum

W65-14

Barraquer


Barraquer lid speculum, .035" Stainless Wire, This design is intended to lay over the temporal, 14mm Closed square blades, 18mm Blade spread, 39mm total length, Stainless



W65-14H

Barraquer


Barraquer lid speculum, .039" Heavy stainless wire, This design is intended to lay over the temporal, 14mm Closed square blades, 20mm Blade spread, 40mm total length, Stainless



W67-14C

McIntyre


McIntyre lid speculum, .035" Stainless Wire, This design is intended to lay over the temporal, 14mm Closed V-Shaped Blades, 19mm Blade spread, 50mm total length, Stainless



W67-140

McIntyre


McIntyre lid speculum, .035" Stainless Wire, This design is intended to lay over the temporal, 14mm Open V-Shaped blades, 19mm Blade spread, 50mm total length, Stainless



W69N-14

Barraquer

Kratz-Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the bridge of the nose, 14mm Open round blades, 18mm Blade spread, 41mm total length, Stainless




Wire Speculum

W69N-14H

Barraquer


Nichamin lid speculum, .039" Heavy stainless wire, This design is intended to lay over the bridge of the nose, 14mm Open Round Blades with wings, 18mm blade spread, 41mm total length, Stainless



W80-10F

Kratz-Barraquer


Kratz-Barraquer Pediatric Speculum, .030" Stainless wire, This design is intended to lay over the temporal, 10mm Open Square Blades with flared tips, 19mm Blade Spread, 35mm total length



W80-14

Kratz-Barraquer


Kratz-Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the temporal, 14mm Open square blades, 18mm Blade spread, 41mm total length, Stainless



W80-14F

Barraquer


Kratz-Barraquer lid speculum, .035" Stainless wire, This design is intended to lay over the temporal, 14mm Open square blades with flared tips, 18mm Blade spread, 41mm total length, Stainless



W80-14FH

Barraquer

Kratz-Barraquer lid speculum, .039" Heavy stainless wire, This design is intended to lay over the temporal, 14mm Open square blades with flared tips, 20mm Blade spread, 41mm total length, Stainless




Wire Speculum

W80N-10

Barraquer


Kratz-Barraquer Pediatric lid speculum, .030" Stainless wire, This design is intended to lay over the bridge of the nose, 10mm Open square blades, 19mm Blade spread, 35mm total length, Stainless



W80N-14

Lehner


Lehner Eye Speculum, .035" Stainless wire, This design is intended to lay over the bridge of the nose, 14mm Open square blades, 18mm Blade spread, 41mm total length, Stainless



W85-10

Barraquer


Barraquer Child lid speculum, .035" Stainless Wire, This design is intended to lay over the temporal, 10mm Solid blades, 15mm Blade spread, 32mm total length, Stainless



W85-14

Barraquer


Barraquer Eye Speculum, .039" Heavy stainless wire, This design is intended to lay over the temporal, 14mm Solid Blades Rounded, 18mm Blade Spread, 40mm total length



W85N-14

Barraquer

Barraquer Eye Speculum, .039" Heavy stainless wire, This design is intended to lay over the bridge of the nose, 14mm Solid Blades Rounded, 18mm Blade Spread, 40mm total length





Diamatrix Ltd

210 Nursery Road
The Woodlands, TX 77380 USA

Tel: 281.367.8081

U.S. Toll Free: 800.867.8081

Fax: 281.292.5481

cs@diamatrix.com

diamatrix.com

